HYDRO-SURG™ Plus Laparoscopic Irrigator with NEZHAT-DORSEY™ SMOKEVAC™ Trumpet Valve and Cojoined Suction and Irrigation Tubing

Instructions for Use

Do not resterilize

Do Not Re-use

STERILE EO

R_x Only

Caution: Read all instructions prior to use.

WARNING: After use, this product is a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws (U.S.A.) and regulations.

CONTRAINDICATIONS:

Use of this device for hysteroscopy or for cavity distention is contraindicated.

Use of this device for intraabdominal irrigation is contraindicated whenever laparoscopy is contraindicated. See operator's manual of your laparoscope for absolute and relative contraindications.

Use of electrosurgery instruments and electrode inserts with an uninsulated metallic tip is contraindicated.

Reorder Numbers:

REF 0026870

HYDRO-SURG™ Plus Laparoscopic Irrigator with NEZHAT-DORSEY™ SMOKEVAC™ Trumpet Valve, with 33cm Tip and Cojoined Suction and Irrigation Tubing

REF 0026880

HYDRO-SURG™ Plus Laparoscopic Irrigator with NEZHAT-DORSEY™ SMOKEVAC™ Trumpet Valve, Handle only, with Cojoined Suction and Irrigation Tubing

MEDICAL EQUIPMENT WITH RESPECT TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1 (2003), ANSI/AAMI ES60601-1 (2005), AND CAN/CSA-C22.2 NO. 60601-1 (2008)
Product contains no user serviceable parts.

BAIRD

DAVOL INC.

PK3797847
143R

Intended Use:

The HYDRO-SURG™ Plus Irrigator is designed to be used with the DAVOL™ Trumpet Valve and tip to provide controlled powered irrigation to, and aspiration of fluids/smoke from, the operative site during laparoscopic surgical procedures (e.g., laparoscopic cholecystectomy and laparoscopic gynecological procedures). It may also be used for resection of filmy adhesions (i.e., hydrodissection) and peritoneal lavage. Appropriate fluids include those which have a specific gravity of approximately 1 (e.g. saline).

Additional specialty probe tips are available from DAVOL, including insulated and non-conductive tips for use with electrosurgery.

Suction/irrigation probe tips, when included, have a non-reflective surface to help avoid laser beam reflection.

With proper adapter, a laser fiber and laser fiber sheath can be passed through the probe handle and tip. Pathology and/or surgeon's choice will dictate preference of suction/irrigation probe tip configuration.

Product Description:

The HYDRO-SURG™ Plus motor driven impeller pump is powered by eight alkaline batteries, size AA, 1.5 volts each, which have an expected shelf life of four years. The batteries, motor, and impeller are all located inside the pump casing.

Contraindications:

- Use of this device for hysteroscopy or for cavity distention is contraindicated.
- Use of this device for intra-abdominal irrigation is contraindicated whenever laparoscopy is contraindicated. See operator's manual of your laparoscope for absolute and relative contraindications.
- Use of electrosurgery instruments and electrode inserts with an uninsulated metallic tip is contraindicated.

Warnings:

- **This device has been designed for single use only. Reuse, reprocessing, resterilization or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device and may lead to device failure which may result in injury to the patient.**
Reuse, reprocessing, resterilization or repackaging may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient or end user.
- **Do not resterilize.**
- **Do not incinerate the batteries as this can cause a risk to the environment.**
- **Only remove batteries for disposal when procedure is completed. Batteries should never be replaced, recharged, put in backwards, or mixed with used or other battery types as these actions may cause explosion or leakage leading to personal injury or can cause a risk to the environment.**
- **As with any surgical procedure where extensive irrigation is utilized, unintended fluid uptake levels should be closely monitored by the attending physician.**
- **When using electrosurgery tips or electrode inserts, follow all instructions, warnings, precautions, and contraindications provided with these devices.**
- **Do not use in an oxygen-enriched atmosphere due to explosion hazard.**
- **Do not occlude suction tip side holes. If side holes become occluded, negative pressure at the tip may increase and cause tissue to be pulled into the tip potentially damaging tissue.**
- **If HYDRO-SURG™ Plus Irrigator is used stacked with or adjacent to other equipment, HYDRO-SURG™ Plus Irrigator should be observed to verify normal operation.**

Precautions:


- If package is damaged or open, do not use product.
- The pump motor runs continuously when fluid is **NOT** present in the pump housing. **Prolonged operation without fluid in the system may damage the pump motor and must be avoided.** This could occur in two ways: 1.) if the pump is turned on prior to priming, turn pump "OFF" and see Step 5 for Instructions for Gravity Priming the Pump; or, 2.) if air enters the pump after bag runs out of fluid, see Step 7 for Instructions for Repriming System in Case of Bag Change.
- Attach pump to I.V. pole before spiking irrigation bag since spike will not securely support pump while in use.
- Use of bottled irrigant is **NOT** recommended. Although bottles can be used, performance is significantly diminished by air that is pulled into the pump before it can be vented from the bottle.
- Adhesions, anatomical anomalies or other conditions may prevent clear visualization of internal anatomy. Do not perform laparoscopic procedures if the internal anatomy cannot be positively identified.
- Laparoscopic procedures should be performed only by persons having adequate training and familiarity with laparoscopic techniques. Consult the medical literature relative to techniques, complications and hazards prior to performance of any laparoscopic procedure.
- This product is for single use and is provided sterile. **DO NOT REUSE.**
- Follow all of the recommended electrosurgical and laser manufacturers' instructions and warnings while operating in these modes.
- Install medical equipment according to the EMC information provided.
- Mobile RF communications equipment can effect medical equipment.

Warnings and Precautions:**WHEN PASSING LASER FIBERS THROUGH PROBE:**

- Do not use CO₂ laser energy systems as they are not compatible for use through the probe handle and tip.
- A thorough understanding of the principles and techniques of laser surgical procedures is necessary to avoid injury to the patient, operator and operating room personnel.
- Pathology and the surgeon's choice will dictate laser tip configuration.
- Proper eye wear is mandatory before activating the laser.
- Use of this device with laser equipment is indicated only in those procedures for which the selected laser fiber and laser system being used have been approved by the FDA.
- To avoid damage to the laser fiber, do not expose tip until suction/irrigation probe has passed through the laparoscopic trocar.
- To avoid injury to the patient do not activate laser until tip is fully inserted and under visual control.


Storage:

- Ambient Temperature: 15°C to 30°C (approx. 59°F to 86°F)
- Relative Humidity: 30% to 65% (Non-Condensing)

**Instructions For Use:**

1. Aseptically remove probe with tubing set from package and inspect for obvious damage. Do not use if damage is evident.

Figure 1

Additional Equipment Required:

- Irrigation Bags
- Suction Irrigation Probe Tip (supplied with REF 0026870)

2. Attach HYDRO-SURG™ Plus Irrigator to I.V. pole.

Note:

- Pump must be attached to the I.V. pole since spike will not securely support pump while in use.
- For non-standard I.V. poles, an optional I.V. Pole Adapter (REF 0026810) is available. Adapter fits 3/8" to 1" (10 mm to 25 mm) diameter IV-poles.

3. Remove spike cap from irrigation spike and connect Irrigator Pump spike to irrigation bag, ensuring pump is at least one foot above patient height. (ref fig.2)

Note:

- Optional Dual Spike Adapter (REF 0026110) is available..

4. Connect suction tubing to suction source. Peel away from irrigation tubing as needed.
5. Check and tighten the knurled Quick Disconnect Adapter at the front of the Trumpet Valve handpiece.
6. Attach selected probe tip to trumpet valve. If using REF 0026870, remove protective sleeve from disposable probe tip prior to use. Ensure that the probe tip is fully seated (on some tips a double click can be felt and/or heard). Additional specialty tips are available from DAVOL, including insulated and non-conductive tips for use with monopolar electrosurgery inserts. Probe tips are quickly and easily exchanged at any time. Always inspect probe tips for damage. Do not use if damage is evident.
7. To prime, depress irrigation button (colored button) and gravity prime unit until all air is out of the system and only irrigant exits the probe. Prime system with probe outside of patient. Rotate the power button to the "I" or "ON" position.
8. Insert shaft of probe through a laparoscopic cannula advancing gently to penetrate the trocar valve. Under appropriate visualization, advance until tip is at desired location.
9. Depress irrigation button (colored button) to irrigate. Depress the unmarked or white button on the Trumpet Valve to suction.

Note: TO REPRIME SYSTEM IN CASE OF BAG CHANGE: If the irrigation bag inadvertently runs out of fluid, system may be reprimed by rotating the power button to the "O" or "OFF" position, spiking a new irrigation bag and repeating steps 5 through 7.
10. **TO ACTIVATE SMOKEVAC™ FEATURE:** Rotate SMOKEVAC™ mechanism counterclockwise to achieve desired rate of smoke/plume removal. SMOKEVAC™ feature will remain at activated setting even after suction valve is depressed and released.

Note: An immediate medium-to-fast setting can be achieved by rotating the SMOKEVAC™ mechanism counterclockwise until it stops and then fully depressing the suction button.
11. **TO DEACTIVATE SMOKEVAC™ FEATURE:** Rotate mechanism clockwise until a complete stop is reached and all incremental markings are visible.
12. Always turn pump "OFF" and remove batteries prior to disposal.

After Use:


1. Ensure Suction tubing is empty and then disconnect Suction connection from Suction Source.
2. Ensure irrigation tubing is empty and then disconnect spike from irrigation bag.

Note: Failure to remove fluids from all tubing may lead to spillage.

Battery Removal:

The batteries, which power the product, may be disposed of separately by pushing in the tab located on the front of the battery casing, pulling the battery casing downward, and removing the batteries from the casing. If the battery casing does not separate easily, carefully insert a sturdy instrument between the pump top and battery casing to aid in prying the pump housing apart.

Figure 2

EMC TABLES

Guidance and manufacturer's declaration – electromagnetic emissions

The HYDRO-SURG™ Plus is intended for use in the electromagnetic environment specified below. The customer or the user of the HYDRO-SURG™ Plus should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The HYDRO-SURG™ Plus uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The HYDRO-SURG™ Plus is suitable for use in all establishments other than domestic.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacturer's declaration – electromagnetic immunity

The HYDRO-SURG™ Plus is intended for use in the electromagnetic environment specified below. The customer or the user of the HYDRO-SURG™ Plus should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	N/A	N/A	Not Applicable - Battery powered
Surge IEC 61000-4-5	N/A	N/A	Not Applicable - Battery powered
Voltage dip, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	N/A	N/A	Not Applicable - Battery powered
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration – electromagnetic immunity

The HYDRO-SURG™ Plus is intended for use in the electromagnetic environment specified below. The customer or the user of the HYDRO-SURG™ Plus should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	N/A 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the HYDRO-SURG™ Plus, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance Battery Operated $d = 1.17 \sqrt{P}$ $d = 1.17 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.33 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- ^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HYDRO-SURG™ Plus is used exceeds the applicable RF compliance level above, the HYDRO-SURG™ Plus should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the HYDRO-SURG™ Plus.
- ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the HYDRO-SURG™ Plus

The HYDRO-SURG™ Plus is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the HYDRO-SURG™ Plus can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the HYDRO-SURG™ Plus as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P}$	80 MHz to 800 MHz $d = 1.17 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.33 \sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.370	0.370	0.737
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.7	11.7	23.3

For transmitters rated at maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

SAFETY INFORMATION AND SYMBOLS DEFINITION

	Danger: Explosion Hazard. Do not use in the presence of flammable anesthetics.	IPX1	Protected against ingress of dripping water.
	Do Not Re-use		Do not resterilize
	Consult Instructions for Use		Type B applied part
	Contents		Contains phthalates
	Do not use if package is damaged	R_X Only	U. S. Federal law restricts this device to sale by or on the order of a physician.

 Manufacturer:
Davol Inc.
Subsidiary of C. R. Bard, Inc.
100 Crossings Boulevard
Warwick, RI 02886 USA
1-401-825-8300 • 1-800-556-6275
Assembled in Mexico

Bard, Davol, Hydro-Surg, Nezhat-Dorsey and SmokEvac are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 2003, 2008, 2010, 2011, 2014 C. R. Bard, Inc. All Rights Reserved.

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Copyright © 2003, 2008, 2010, 2011, 2014 C. R. Bard, Inc. Tous droits réservés.
Daval Inc.