Reorder Number

REF **5551000**

Non-sterile

Caution: Read all instructions prior to use.

Rx Only

X **X-STREAM™**
Laparoscopic Irrigation Controller

Operator's Manual

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Note: The X-STREAM™ Laparoscopic Irrigation Controller will show no visual or audible signs of powering on until the X-STREAM™ Laparoscopic Irrigation Tubing Set is loaded into the Controller (even when the Controller is plugged into an electrical receptacle).

DAVOL INC.

PK3790414
Rev. 2018/01

Released

	Safety Information and Symbols Definition	en 1
	Accessories	en 1
Section 1	Introduction	en 2
Section 2	Warnings	en 2
	A. X-STREAM™ Laparoscopic Irrigation Controller	
	B. X-STREAM™ Laparoscopic Irrigation Tubing Set	
Section 3	Precautions	en 2
Section 4	Specifications	en 3
Section 5	User Interface (Description/Function)	en 3
Section 6	Controller Setup and Operating Instructions	en 4-5
Section 7	Cleaning and Maintenance	en 5
Section 8	Calibration/Verification Testing	en 5-6
	A. Grounding Resistance	
	B. System Calibration Procedure	
	C. Conclusions	
Section 9	Troubleshooting	en 6
Section 10	Controlling System Audio Indicator	en 7
Section 11	Appendix	en 8
	A. Return Policy Notification	
	B. Limited Product Warranty; Disclaimer and Limitation of Liability	
	C. Product Changes	
	D. Force Majeure	
	E. Safety Agency Testing	
	Product Illustration	en 9

SAFETY INFORMATION and SYMBOLS DEFINITION

**CAUTION: Read all instructions prior to use.**

The exclamation point within an equilateral triangle is intended to alert the user to the presence of important operating and maintenance (servicing) instructions in the literature accompanying the device.

The lightning flash with arrowhead symbol, within an equilateral triangle, is intended to alert the user of the presence of uninsulated "dangerous voltage" within the product's enclosure that may be of sufficient magnitude to constitute a risk of electric shock to persons.

Maximum flow rate setting

Other Patents Pending

Type B Applied Part

Default flow rate setting

Fuse

IPX4

Protected against the effects of splashing water.

PRIME Prime

Reusable

Warranty void if seal is broken.

"WARNING": Risk of fire.
Replace fuse as marked.

Contents

U.S. Federal law restricts this device to sale by or on the order of a physician.

Indicates separate collection for electrical and electronic equipment.

Biological Risks

Temperature Limitation

Consult Instructions for Use

Protective earth (ground)

Non-Sterile

Contains phthalates

ACCESSORIES

REF 5552000

X-STREAM™ Laparoscopic Irrigation Tubing Set with SMOKEVAC™ Trumpet Valve, 5mm x 33cm Probe Tip with Holes and Cojoined Suction/Irrigation Tubing

REF 5552002

X-STREAM™ Laparoscopic Irrigation Tubing Set with NEZHAT-DORSEY™ Trumpet Valve, 5mm x 33cm Probe Tip with Holes and Cojoined Suction/Irrigation Tubing

REF 5552001

X-STREAM™ Laparoscopic Irrigation Tubing Set with SMOKEVAC™ Trumpet Valve and Cojoined Suction/Irrigation Tubing, without Probe Tip

REF 5552003

X-STREAM™ Laparoscopic Irrigation Tubing Set with NEZHAT-DORSEY™ Trumpet Valve and Cojoined Suction/Irrigation Tubing, without Probe Tip

The X-STREAM™ Laparoscopic Irrigation Controller is a versatile and reusable electromechanical controller (Figure 1, Page en 3). It is small, lightweight, and pole-mounted. It has two easily adjustable flow rate settings, and may be used with any of the DAVOL™ Disposable X-STREAM™ Laparoscopic Irrigation Tubing Sets.

SECTION 2 WARNINGS

A. X-STREAM™ Laparoscopic Irrigation Controller

1. The performance of the X-STREAM™ Laparoscopic Irrigation Controller depends on proper setup. Please refer to tubing set Instructions for Use for setup instructions.
2. Use of any modified or unauthorized disposables with the X-STREAM™ Laparoscopic Irrigation Controller is contraindicated and may result in undesirable system performance.
3. If performing electrical safety testing follow IEC 60601-1:2006 recommendations.
4.  The X-STREAM™ Laparoscopic Irrigation Controller contains internal components that could cause electrical shock hazard. Do not open the housing of the Controller or attempt to service the system.
5. The power supply cord is not rewirable. Contact DAVOL™ Customer Service if you experience any functional issues with the X-STREAM™ Laparoscopic Irrigation Controller.
6. To avoid the possibility of electrical shock, always remove the tubing set and disconnect the power cord prior to cleaning.
7. Do not use the X-STREAM™ Laparoscopic Irrigation Controller in the presence of a FLAMMABLE ANESTHETIC MIXTURE WITH AIR or with OXYGEN OR NITROUS OXIDE.
8. Do not autoclave or steam sterilize the X-STREAM™ Laparoscopic Irrigation Controller or subject the pump to temperatures in excess of 60°C. To avoid possibility of electric shock, always disconnect the power cord prior to cleaning.
9. Do not touch the X-STREAM™ Laparoscopic Irrigation Controller and the patient simultaneously.
10. No user serviceable items. Contact DAVOL™ Customer Service if you experience any functional issues with the X-STREAM™ Laparoscopic Irrigation Controller.
11. The tubing set will provide gravity flow if the pumping chamber is not properly loaded into the controller's pumping chamber receptacle.
12. Do not modify this equipment as it may compromise the safety and effectiveness of the device.
13. If the X-STREAM™ Laparoscopic Irrigation Controller is installed in the vicinity of a RF communication device (see tables) proper operation of the X-STREAM™ Laparoscopic Irrigation Controller should be observed.
14. The X-STREAM™ Laparoscopic Irrigation Controller should be installed in accordance with the EMC information provided in this document.
15. Do not use accessories other than those declared by DAVOL™ as this may negatively impact X-STREAM™ Laparoscopic Irrigation Controller's emissions and immunity.

B. X-STREAM™ Laparoscopic Irrigation Tubing Set

1. As with any surgical procedure where extensive irrigation is utilized, intravasation levels should be closely monitored by the attending physician.
2. The pressure accuracy of the X-STREAM™ Irrigation System depends on proper setup. It is important that the bottom of the X-STREAM™ Laparoscopic Irrigation Controller be positioned at least 12 inches above the highest height of the patient to provide sufficient pressure for optimal performance of the system.
3. The X-STREAM™ Irrigation System is capable of generating high irrigation impact forces and flow rates. Use the appropriate, required pressure necessary at all times to accomplish the procedure.
4. The X-STREAM™ Irrigation System will generate even higher irrigation impact forces when used with tips that have an irrigation channel smaller than 5mm (compared to irrigation impact forces generated with a 5mm standard probe tip).
5. Electrosurgery instruments and electrode inserts must not be used with uninsulated metallic probe tips. Electrosurgery instruments and inserts should be used only with available DAVOL™ non-conductive or insulated probe tips.
6. When using electrosurgical attachments or inserts, follow all instructions, warnings, precautions, and contraindications provided with these devices.

SECTION 3 PRECAUTIONS

1. Read all the instructions provided with the X-STREAM™ Laparoscopic Irrigation Tubing Set prior to use.
2. The X-STREAM™ Laparoscopic Irrigation Controller automatically turns "ON" when the tubing set is loaded into the controller. The "DAVOL" logo and "PRIME" LED's will light to indicate that the controller is on. There is no "ON/OFF" switch.
3. Do not remove the square, metal contact label located at the top back of the tubing set pumping chamber.
4. Solvents should not be used on the controller's logo and keypad as they could cause damage to the device.
5. When air enters the system, the controller returns to the "PRIME" state, the "PRIME" text will blink, and the audio indicator will sound.
6. Gravity prime and perform all priming with the probe tip outside of patient.
7. Please use special precautions regarding the use of the X-STREAM™ Irrigation System. The X-STREAM™ Irrigation System is to be used in EMC environments that are suitable for Group 1, Class A medical equipment.
8. Grounding reliability can only be achieved when the power cord supply is connected to an equivalent receptacle marked "Hospital Only" or "Hospital Grade".
9. The X-STREAM™ Laparoscopic Irrigation Controller should be positioned so that it is not difficult to disconnect device from the power source, in order to facilitate access.
10. The X-STREAM™ Laparoscopic Irrigation Controller should not be placed adjacent to or stacked with other equipment unless nominal operation of the X-STREAM™ Irrigation System has been observed as it may cause a malfunction of the device.

Physical	Electrical	System	Environmental/Storage
Width: 16.0cm	Voltage: 100-240VAC	Classification: Class 1	Generated Noise: 85dBA
Length: 15.2cm	Frequency: 50/60 Hz	Mode: Intermittent Operation	Ambient Temperature: 15°C to 30°C
Height: 18.4cm	Fuse(s): T 2,5AH, 250V	Type:  Type B Applied Part	Relative Humidity: 30% - 65% (Non-Condensing)
Weight: 1.2kg		Ingress Protection: IPX4	Atmospheric Pressure: 70kPa to 100kPa
		Continuous Irrigation Time (max): 15 sec/cycle	

SECTION 5 USER INTERFACE

ITEM	DESCRIPTION	FUNCTION
1	“DAVOL” logo LED	An illuminated logo indicates the presence of a properly loaded disposable tubing set and that the controller is “ON”, waiting for “PRIME” or “READY”.
2	“PRIME” LED	Illuminated “PRIME” text indicates that the tubing set needs to be primed.
3	“HI” setting LED button	Depressing this button activates the maximum flow rate. There are only two flow rate settings.
4	“LO” setting LED button	Depressing this button activates the default flow rate. There are only two flow rate settings.
5	IV pole clamp	Mounts controller to standard IV poles.
6	Fuse holder	Holds two T 2,5AH, 250V fuses. Lift tab on inside of IEC outlet to access fuse holder.
7	IEC power cord	Attachable power cord compatible with IEC power cord receptacle. Attaches to (8).
8	IEC power cord receptacle	Female receptacle for IEC compatible, attachable power cords.
9	Power cord retention strap stud	For securing the strap.
10	Power cord retention strap	Useful in retaining the power cord during storage of the controller.
11	Metal case	 Grounded surface for testing leakage current.
12	Pumping chamber receptacle	Receptacle to receive disposable tubing set pumping chamber.

Figure 1



Figure 2

Important: Thoroughly review all of the operating instructions in this manual and in the X-STREAM™ Laparoscopic Irrigation Tubing Set Instructions for Use before using the X-STREAM™ Laparoscopic Irrigation Controller.

X-STREAM™ Laparoscopic Irrigation Controller Setup

1. With the bottom of the X-STREAM™ Laparoscopic Irrigation Controller at least 12 inches above the highest height of the patient at any time during surgery, mount the controller on the IV pole. To mount, loosen pole clamp knob by turning it. Position the pole clamp around the IV pole and tighten the knob to secure the controller to the IV pole.
2. Ensure the X-STREAM™ IEC/Hospital Grade power cord is plugged into the IEC receptacle on the back of the controller. Also ensure that the power cord is connected to an operable Hospital Grade AC outlet.
3. Hang irrigation bag(s) on the IV pole. Because of the high flow rates generated by the X-STREAM™ Irrigation System, 2000cc or larger IV bag sizes are recommended to minimize bag changes.
4. If spiking two bags is desirable, use the optional Dual Spike Adapter (REF 0026110).

Note: Changing the setting is possible both with and without the motor running.

Note: The flow rate setting cannot be changed when the controller is in the “PRIME” mode.

Precaution: Do not remove the square, metal contact label located at the top back of the disposable tubing set pumping chamber.

Opening the Disposable Tubing Set Package

1. The circulating nurse opens the package and delivers the tubing set onto the sterile field.
2. The scrub nurse retains the trumpet valve handpiece and checks and tightens the knurled Quick Disconnect Adapter at the front and the knurled knob plug at the back of the Trumpet Valve.
3. The scrub nurse securely attaches a probe tip by pushing the probe tip over the Quick Disconnect Adapter (on some tips a double “click” is heard), then hands the pumping chamber (X-STREAM™ clear cartridge) to the circulating nurse. Refer to the X-STREAM™ Laparoscopic Irrigation Tubing Set Instructions for Use for information on the probe handle and tips.
4. The circulating nurse connects the suction tubing to the suction source. Peel away the suction tubing from the irrigation tubing as needed.

Loading the Disposable Tubing Set

Precaution: Read all directions for the X-STREAM™ Laparoscopic Irrigation Tubing Set before beginning the loading procedure.

1. Slide the pumping chamber downward into the controller until it snaps into place flush with the top and face of the controller.
Precaution: The controller automatically turns “ON” when the tubing set is loaded into the controller. The “Davi” logo and “PRIME” LED’s will light to indicate that the controller is “ON”. There is no “ON/OFF” switch.
2. Remove spike cap and spike the irrigation bag. Leave sufficient slack in line to facilitate bag placement/removal. If spiking two bags is desirable, use the optional Dual Spike Adapter (REF 0026110).

Priming the System

Note: The system includes audible and visual indicators to let the user know priming is required. The “PRIME” LED light will flash and the audible indicator will sound when the system detects the presence of the pumping chamber upon initial setup. See Section 10, page en 7, for controlling the audio indicator.

Precaution: Gravity prime and perform all priming with probe tip outside of patient.

1. To prime, depress irrigation button (colored button) until all air is out of the system and only irrigant exits the probe tip.
2. Release the irrigation button.
3. If the pumping chamber is not primed properly, the “PRIME” text will continue to flash and the audio indicator will continue to sound after the irrigation button is released.

Note: Removal and reloading of the pumping chamber will return the controller to the “PRIME” mode. The “DAVOL” logo LED will illuminate to indicate that the X-STREAM™ Laparoscopic Irrigation Controller is “ON”. The “PRIME” text will flash and the audio indicator will sound to indicate that the tubing set needs to be primed.

4. When ready to irrigate, depress the irrigation button (colored button) on the trumpet valve handpiece.

Note: The motor will run when the irrigation button is depressed and will stop after it is released. A safety feature turns the pump off after 15 seconds of continuous use. Release the irrigation button and then re-depress to reactivate the pump.

Basic Operation

1. Changing the Flow Rate Setting

1. There are two flow rate settings available on the X-STREAM™ Laparoscopic Irrigation Controller. The LED lights to identify the activated flow rate setting.
2. The controller automatically activates at the low setting. To select the maximum flow rate setting, the circulating nurse should press the “HI” button.
3. Switch back and forth between the “LO” and “HI” settings throughout the procedure by pressing the appropriate button.
4. The surgeon can control the flow at the handpiece.

**Default
Flow Rate**
2500 - 3800 cc/min

**Maximum
Flow Rate**
3400 - 5000 cc/min

Note: Changing the setting is possible both with and without the motor running.

Note: The flow rate setting cannot be changed when the controller is in the “PRIME” mode.

Note: The NEZHAT-DORSEY™ Trumpet Valves offer variable flow control. For example, depressing the irrigation button on the handpiece halfway delivers approximately one-half of the available irrigation pressure/flow.

2. Changing Irrigation Bags

Note: The system includes audible and visual indicators to let the user know priming is required. The “PRIME” and “HI” LED lights will blink and the audio indicator will sound when the bag runs out of fluid.

Precaution: Perform all priming with the probe tip outside the patient. If the bag(s) runs out of irrigant, the motor will stop running. When air enters the system, the controller returns to the “PRIME” state, the “PRIME” and “HI” text will blink, and the audio indicator will sound.

1. Hang a new bag(s) and spike.
2. Depress the irrigation button (colored button) on the handpiece until fluid fills the pumping chamber, and then, while continuing to hold the irrigation button down, press the “HI” button on the controller to power prime the system. Release irrigation button after 3 seconds.

Note: The “PRIME” text and the audio indicator will go off and the flow rate setting LED will light up. The system returns to the flow rate setting in use prior to bag runout.

If the “PRIME” text remains lit, and the audio indicator is sounding, the tubing set has not been primed properly. Continue power priming to remove all air from the system.

3. Removing the Disposable Tubing Set

1. If using a Dual Spike Adapter (REF 0026110) close the appropriate clamp(s) on the bag spike arms of the disposable tubing set.
2. Remove suction tube from suction device.
3. Remove the pumping chamber from the controller by sliding it upwards.
4. Discard the disposable tubing set in accordance with hospital waste disposal practices.

Precaution: After use, the tubing set, trumpet valve handpiece and probe tip are a potential biohazard. Handle and dispose of as required by hospital policy and applicable laws.

5. Clean the X-STREAM™ Laparoscopic Irrigation Controller and the X-STREAM™ Power Cord after each use (See **Section 7, Cleaning and Maintenance**).

SECTION 7 CLEANING AND MAINTENANCE

WARNING: Do not autoclave or steam sterilize the X-STREAM™ Laparoscopic Irrigation Controller or power cord. Do not subject the controller or power cord to temperatures in excess of 60°C. To avoid possibility of electric shock, always disconnect the power cord from the power supply outlet prior to cleaning and maintenance.

1. Cleaning

Note: Clean the X-STREAM™ Laparoscopic Irrigation System Controller and Power Cord after each use. Follow the pre-saturated disposable wipe manufacturer’s instructions for use.

- A. Cleaning the X-STREAM™ Laparoscopic Irrigation System Controller:
 1. Use a pre-saturated germicidal disposable wipe (i.e. PDI Super SANI-CLOTH® Germicidal Disposable Wipe) that is intended for cleaning and disinfection in health care environments to thoroughly wipe the Surface of the X-STREAM™ Controller (including the pumping chamber receptacle, see **Section 5, USER INTERFACE, Item 12**).
 2. Repeat Step 1 with a new wipe, ensuring all visible soil and debris is removed.
 3. If any visible soil or debris remains on the device, use as many additional wipes as necessary to ensure that it is completely removed.
- B. Cleaning the X-STREAM™ Laparoscopic Irrigation System Controller Power Cord:
 1. Use a fresh pre-saturated germicidal disposable wipe (i.e. PDI Super SANI-CLOTH® Germicidal Disposable Wipe) that is intended for cleaning and disinfection in health care environments to thoroughly wipe the X-STREAM™ Laparoscopic Irrigation Controller Power Cord.
 2. If any visible soil or debris remains on the device, use as many additional wipes as necessary to ensure that it is completely removed.
- C. Store the X-STREAM™ Laparoscopic Irrigation Controller and Power Cord for future use.

2. Maintenance

Fuse Replacement

- A. In the event that fuses need replacement, perform the following steps:
 1. Locate the fuse holder, (See **Section 5, USER INTERFACE, Item 6** for accessing fuses) lift tab on inside of IEC outlet to access fuse holder
 2. Remove both fuses and replace with two new T 2,5AH, 250V fuses.
 3. Reinsert the fuse holder with new fuses into the IEC outlet

SECTION 8 CALIBRATION / VERIFICATION TESTING

A. Grounding Resistance

Power cord integrity may be checked without opening the controller casing by measuring grounding resistance from the Power Entry Module Metal Casing on the back of the controller to the ground pin of the attachment plug per NFPA 99 Standard for Health Care Facilities (1996 edition) Chapter 7-5.1.3.2 Resistance. All other surfaces are isolated from ground. Perform electrical safety testing following IEC 60601-1:2006 recommendations.

B. System Calibration Procedure

Electrical Voltage Verification Protocol


Purpose

To verify that the reusable controller is delivering appropriate voltage to the disposable tubing set.

Equipment
RMS Multimeter

Procedure

1. Set up disposable tubing set and prime system according to Instructions for Use.
2. Maintain a full 3L IV bag at the height described in the Instructions for Use for the disposable tubing set.
3. On the RMS Multimeter, connect conductor leads to Common (-) and Volts/Ohms (+).
4. Set dial on RMS Multimeter to measure DC Volts (\bar{V})
5. Contact the probed end of the conductor leads to the metal contacts at the bottom of the pumping chamber section of the disposable tubing set as shown in the **Electrical Test Setup Diagram**.
6. Activate the **(LO)** button on the controller. Irrigate and record maximum voltage reading on RMS Multimeter while irrigating.
7. Activate the **(HI)** button on the controller. Irrigate and record maximum voltage reading on RMS Multimeter while irrigating.

Electrical Test Setup Diagram

Calibration Table

Controller Setting	Acceptable Range of DC Voltage (± 1VDC)
DEFAULT	19 VDC
MAXIMUM	23 VDC

C. Conclusions

Any deviations from the test results presented above should be reported to the customer service department of the manufacturer, Davol Inc., 100 Crossings Boulevard, Warwick, RI 02886, 1-800-556-6275. If outside of the U.S.A., deviations should be reported to your local X-STREAM™ System distributor.

SECTION 9 TROUBLESHOOTING

PROBLEM	POSSIBLE SOLUTIONS / ACTIONS
"DAVOL" Logo LED's do not illuminate when the disposable tubing set pumping chamber is installed into pumping chamber receptacle of the controller.	<ul style="list-style-type: none"> • Connect power cord securely to the IEC power cord receptacle on the back of the controller. • Connect power cord securely to an operable Hospital Grade AC outlet. • Load the disposable tubing set pumping chamber into the pumping chamber receptacle of the controller (Press firmly downward until the disposable is flush with the top surface of the controller.) • Check fuses. If replacement is needed, see Section 7, Cleaning and Maintenance for fuse replacement. • Check the back of the disposable tubing set pumping chamber to ensure that the metallic label is installed. • Install new disposable tubing set. • Call DAVOL™ Customer Service if problem persists.
"DAVOL" Logo LED's blink after the disposable tubing set pumping chamber is installed but pump will not start.	<ul style="list-style-type: none"> • System fault detected. Electrical continuity to the disposable pumping chamber may be interrupted. • Install new disposable tubing set. • Clean metal contacts on the pumping chamber receptacle of the controller (See Section 8, CALIBRATION/VERIFICATION TESTING, Electrical Test Setup Diagram for location of contacts). • Call DAVOL™ Customer Service if problem persists.
"PRIME" text blinks after system is gravity primed.	<ul style="list-style-type: none"> • Depress irrigation button on trumpet valve handpiece. Prime unit to remove all trapped air. • Call DAVOL™ Customer Service if problem persists.
"DAVOL" Logo LED's are illuminated but fluid is not flowing through the system.	<ul style="list-style-type: none"> • Check that the irrigation bag(s) are not empty. • Check that the irrigation bag(s) are spiked completely. • Open the clamp(s) for the appropriate irrigation bags. • Remove any kinks in the flexible tubing. Ensure that the fluid path is open and fluid is flowing toward the trumpet valve handpiece. • Check that the appropriate button (colored button) is depressed on the trumpet valve handpiece. • Call DAVOL™ Customer Service if problem persists.
Low pressure and low flow.	<ul style="list-style-type: none"> • Ensure that the disposable tubing set pumping chamber is properly seated into the pumping chamber receptacle of the controller. • Ensure that spike is inserted fully into the irrigation container. • Completely open the clamps on the irrigation bag spike arms. • Remove any kinks in the flexible tubing. Ensure that the fluid path is open and fluid is flowing toward the trumpet valve handpiece. • Install new disposable tubing set. • Call DAVOL™ Customer Service if problem persists.
Irrigation flow does not increase or decrease when flow settings are selected from the keypad.	<ul style="list-style-type: none"> • Ensure that the disposable tubing set pumping chamber is properly seated into the pumping chamber receptacle of the controller. • Completely open the clamps on the irrigation bag spike arms. • Remove any kinks in the flexible tubing. Ensure that the fluid path is open and fluid is flowing toward the trumpet valve handpiece. • Firmly depress irrigation button on the trumpet valve handpiece. • Call DAVOL™ Customer Service if problem persists.
Motor does not run. "DAVOL" Logo and flow setting buttons ("HI" or "LO") are illuminated.	<ul style="list-style-type: none"> • Check that the irrigation bag(s) are not empty. • Remove any kinks in the flexible tubing. Ensure that the fluid path is open and fluid is flowing toward the trumpet valve handpiece. • Ensure that the irrigation button on the trumpet valve handpiece is depressed. • Install new disposable tubing set. • Call DAVOL™ Customer Service if problem persists.
Motor runs intermittently without cause.	<ul style="list-style-type: none"> • Depress irrigation button on the trumpet valve handpiece. Prime unit to remove all trapped air. • Avoid vigorous shaking of the irrigation tubing during use. • Call DAVOL™ Customer Service if problem persists.
Motor runs continuously without cause.	<ul style="list-style-type: none"> • Depress irrigation button on the trumpet valve handpiece. Prime unit to remove all trapped air. • Check for fluid leaking at the trumpet valve handpiece or irrigation tip. • Check for fluid leaking at the pumping chamber. • Call DAVOL™ Customer Service if problem persists.

Please refer to illustration.

A. Silencing System Audio Indicator


- To silence the system audio indicator, press and hold down the “LO” button for 3 seconds, until a single beep sounds.
- Release button after the beep and system audio indicator will turn off automatically.
- Wait 3 seconds before resuming normal button operation.

B. Activating System Audio Indicator

- To activate the system audio indicator, press and hold down the “LO” button for 3 seconds, until a single beep sounds.
- Release button after the beep and system indicator alarm will turn on automatically.
- Wait 3 seconds before resuming normal button operation.

C. Adjusting the Volume of the System Audio Indicator

- To adjust the volume of the system audio indicator, hold down the “LO” button for 3 seconds, until a single beep sounds.
- Release button after beep and press the “HI” button.
- Continue to press button until desired volume is achieved.
- Once at desired volume, release the button. This will automatically set your selection.
- Wait 3 seconds before resuming normal button operation.

EMC TABLES

Guidance and manufacturer's declaration – electromagnetic immunity			
The X-STREAM™ Controller is intended for use in the electromagnetic environment specified below. The customer or the user of the X-STREAM™ Controller should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the X-STREAM™ Controller requires continued operation during power mains interruptions, it is recommended that the X-STREAM™ Controller be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Recommended separation distances between portable and mobile RF communications equipment and the X-STREAM™ Controller			
The X-STREAM™ Controller is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the X-STREAM™ Controller can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the X-STREAM™ Controller as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P}$	80 MHz to 800 MHz $d = 1.17 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33 \sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.370	0.370	0.737
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.7	11.7	23.3
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The X-STREAM™ Controller is intended for use in the electromagnetic environment specified below. The customer or the user of the X-STREAM™ Controller should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the X-STREAM™ Controller, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = 1.17 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.17 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.33 \sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ³ should be less than the compliance level in each frequency range. ³ Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the X-STREAM™ Controller is used exceeds the applicable RF compliance level above, the X-STREAM™ Controller should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the X-STREAM™ Controller.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic emissions		
The X-STREAM™ Controller is intended for use in the electromagnetic environment specified below. The customer or the user of the X-STREAM™ Controller should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The X-STREAM™ Controller uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The X-STREAM™ Controller is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

A. Return Policy Notification

If the X-STREAM™ Laparoscopic Irrigation Controller does not function after troubleshooting (See Section 8) please call Davol Customer Service at 1-800-556-6275 to obtain a Return Material Authorization (RMA) number.

Do not attempt to service the controller. Servicing the controller voids the warranty of the device. The controller contains void seals that will be broken if the device is opened. If the void seal is broken when the controller is returned to Davol, the warranty is void.

Full Credit Returns. Product returns receive full credit (less any rebate paid or credits issued by Davol to customer with respect to such product) provided (i) defective product is returned within applicable warranty period, (ii) the product is subject to a Recall, or (iii) product shipped in error by Davol and received by Davol and is in saleable condition and in original packaging within 30 day of invoice.

Full Credit Returns (All Other). Products returned (other than those outlined above in Full Credit Returns) within 30 business days of invoice date will receive full credit (less any rebate paid or credits issued by Davol to customer with respect to such product). Products returned within 31- 90 days of corresponding Davol invoice date will receive full credit (less any rebate paid or credits issued by Davol to customer with respect to such product) minus a 30% re-stocking fee.

Ineligible Returns. Davol products are ineligible for return if (i) product was shipped in error and returned and received by Davol more than 30 days after the date of the applicable Davol invoice, (ii) Non-defective Davol product returned in a quantity less than Davol's original unit of sale, (iii) Non-defective product returned more than 90 days after the date of corresponding Davol invoice unless part of a Recall (iv) Defective product returned after expiration date of applicable warranty period unless part of a Recall (v) Product sold non-sterile that has been subject to sterilization processing (vi) Product sold for single use that has been re-used or reprocessed (vii) Product that has been altered, further manufactured, packaged, processed, abused, or misused (ix) Product that has been adjusted or repaired by anyone other than an authorized Davol representative, (x) Product that is a "custom" device unless such product is defective for a reason other than manufacture to customer's specifications, or (xi) Products damaged in transit.

RMA. Customer must obtain a Return Material Authorization (RMA) number from Davol or a Davol account representative before any Davol product can be returned for repair, replacement, or credit. Call Davol divisional Customer Service to obtain RMA number. To obtain return RMA number, customer must provide Davol with (i) the Davol catalog number and quantity of Davol product to be returned/repared; (ii) the reason for the return/repair; (iii) a description of the product being returned/repared; (iv) the name and telephone number of a customer contact who may be called if further information is required; (v) at least one of the following: (a) the applicable customer purchase order number (even if product provided at no charge), (b) the corresponding Davol invoice number, and (c) the applicable Davol product lot or serial number. All authorized return products must reference the assigned RMA number and must be shipped freight prepaid by customer with confirmation that the product has been stored according to label storage conditions where applicable. Any authorized returned product sent freight collect will be refused by Davol and returned to customer at their expense. Davol will, at its discretion, refund or credit customer for all freight charges incurred in connection with the return of authorized return products. If product is being returned for repair that is covered by a written limited product warranty, a copy of the original invoice must be sent with the product. The cost of any repair not covered by a written limited warranty must be paid by the customer.

B. Limited Product Warranty; Disclaimer and Limitation of Liability

Davol warrants to the original purchaser that each Davol product set forth in the price list will be free from defects in material and workmanship for the period set forth in the labeling of the particular Davol product or, if no such period is set forth in the labeling, for a period of one (1) year from the date of purchase. Defective Davol product may be returned to Davol for repair, replacement, refund or credit, at Davol's option, in accordance with Davol's return goods and allowance policy. Any alteration, abuse, misuse, further manufacture, packaging, processing, adjustment or repair by any person or entity other than Davol or a person or entity authorized in writing by Davol shall void this limited product warranty ab initio. THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THE LIABILITY AND REMEDY STATED IN THIS LIMITED PRODUCT WARRANTY WILL BE THE SOLE LIABILITY OF Davol, AND REMEDY AVAILABLE TO CUSTOMER, FOR Davol PRODUCTS WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE) OR OTHERWISE, AND Davol WILL NOT BE LIABLE TO CUSTOMER FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE HANDLING OR USE OF DAVOL PRODUCTS EVEN IF DAVOL HAS BEEN ADVISED OF THE POSSIBILITY OR LIKELIHOOD OF SUCH DAMAGES. IN NO EVENT WILL DAVOL BE LIABLE FOR ANY CLAIM, LOSS OR DAMAGE ARISING OUT OF OR RELATING TO, IN WHOLE OR IN PART, ANY PURCHASE ORDER, THESE TERMS AND CONDITIONS OF SALE, OR ANY DAVOL PRODUCT, WHETHER IN CONTRACT, TORT, NEGLIGENCE OR OTHERWISE, IN EXCESS OF THE AMOUNT PAID BY CUSTOMER TO DAVOL PURSUANT TO THE PURCHASE ORDER TO WHICH SUCH CLAIM, LOSS OR DAMAGE RELATES.

C. Product Changes

All products and product specifications identified in the Price List are based upon the information available to Davol at the time of publication. Davol reserves the right to discontinue any product or to change any product specifications without notice.

D. Force Majeure

Davol will not be liable for its failure to perform, or any delay in performance, due to strikes, fire, explosion, flood, riot, lock out, injunction, interruption of transportation, unavoidable accidents, acts of government or a public enemy, terrorism, inability to obtain supplies at reasonable prices, or any other causes beyond its control.


E. Safety Agency Testing

The X-STREAM™ Laparoscopic Irrigation Controller and electrical accessories have been inspected, tested and certified and are in compliance with the applicable requirements of the Standards for Medical Electrical Equipment, IEC 60601-1, UL 60601-1 and CSA C22.2 NO. 601.1.

The X-STREAM™ Laparoscopic Irrigation Controller has been tested and has demonstrated protection against SAFETY HAZARDS caused by harmful ingress of liquids due to splashing, per IEC 60601-1 and UL 60601-1.

MEDICAL EQUIPMENT
WITH RESPECT TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY
IN ACCORDANCE WITH UL 60601-1 (2003), ANSI/AAMI ES60601-1 (2005),
AND CAN/CSA-C22.2 NO. 60601-1 (2008)
Product contains no serviceable parts.

Rx Only

Bard, Davol, SmokEvac, Nezhat-Dorsey, and X-Stream are trademarks and/or registered trademarks of C. R. Bard, Inc. All other trademarks are property of their respective owner(s).

Copyright © 2005, 2008, 2011, 2014, 2018 C. R. Bard, Inc. All Rights Reserved.

Record the Model Number and Serial Number of the X-STREAM™ Laparoscopic Irrigation Controller and retain for future reference.

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Signature Manifest

Document Number: BAW3790009
Title: X-Stream operators manual M/L

Revision: 2

All dates and times are in Eastern Standard Time.

DA CR-12163 Pavel Plugovoy

DC Change Request

Name/Signature	Title	Date	Meaning/Reason
Documentation Placeholder (DOCHOLDER)	Place Holder		
Julia Kue (JKUE)	DOCUMENTATION SPECIALIST		
Les Torrey (LORREY)	ASSOCIATE DIRECTOR, QA		
Brenda Schultz (BSCHULTZ)	DOCUMENTATION CLERK		
Denise Menard (DMENARD)	SUPERVISOR, DOCUMENTATION CTRL	11 May 2018, 11:04:13 AM	Approved

DC Check In 2

Name/Signature	Title	Date	Meaning/Reason
John Golson (JGOLSON)	Pkg Sr. Pkg Engineer		
Pam Fotopoulos (PFOTOPOULOS)	PD - Packaging Coordinator		
Pavel Plugovoy (PPLUGOVOY)	PKG PACKAGING ENG COORDINATOR	14 May 2018, 11:29:25 AM	Complete

Team Review/Approval

Name/Signature	Title	Date	Meaning/Reason
Karen Kullas (KKULLAS)	AME Staff Engineer	14 May 2018, 12:07:10 PM	Approved
Pam Fotopoulos (PFOTOPOULOS)	PD - Packaging Coordinator	14 May 2018, 12:09:03 PM	Approved
John Golson (JGOLSON)	Pkg Sr. Pkg Engineer	15 May 2018, 08:15:45 AM	Approved
Tara Smith (TSMITH)	MARKETING PRODUCT MGR	15 May 2018, 09:51:48 AM	Approved
Pavel Plugovoy (PPLUGOVOY)	PKG PACKAGING ENG COORDINATOR	15 May 2018, 10:56:27 AM	Approved
Gail Dow (GDOW)	RA Sr Reg Affairs Specialist	15 May 2018, 11:09:13 AM	Approved
Christina Liese (CLIESE)	QA AQE I	15 May 2018, 01:52:58 PM	Approved

RA Attachments

Name/Signature	Title	Date	Meaning/Reason
Gail Dow (GDOW)	RA Sr Reg Affairs Specialist	15 May 2018, 02:07:09 PM	Complete

DC Final Review

Name/Signature	Title	Date	Meaning/Reason
Denise Menard (DMENARD)	SUPERVISOR, DOCUMENTATION CTRL	15 May 2018, 02:35:29 PM	Complete
Brenda Schultz (BSCHULTZ)	DOCUMENTATION CLERK	15 May 2018, 02:42:21 PM	Complete

Release

Name/Signature	Title	Date	Meaning/Reason
Denise Menard (DMENARD)	SUPERVISOR, DOCUMENTATION CTRL		
Brenda Schultz (BSCHULTZ)	DOCUMENTATION CLERK	15 May 2018, 02:42:44 PM	Approved

Notify

Name/Signature	Title	Date	Meaning/Reason
Bethanie Benoit (BBENOIT)	MKT Marketing Dir	15 May 2018, 02:42:45 PM	Email Sent
Debra Griffith (DGRIFFITH)	MGR, MEDICAL SERVICES & SUPPOR	15 May 2018, 02:42:45 PM	Email Sent
Pavel Plugovoy (PPLUGOVOY)	PKG PACKAGING ENG COORDINATOR	15 May 2018, 02:42:45 PM	Email Sent

Released

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Copyright © 2005, 2008, 2011, 2014, 2018 C. R. Bard, Inc. Tous droits réservés.
Filiale de C. R. Bard, Inc.