BARD® Mesh Pre-shaped

Monofilament Polypropylene for the
Reconstruction of Tissue Defects

Instructions for Use

Mode d'emploi
Gebrauchsanweisung
Istruzioni per l'uso
Instrucciones de uso
Gebruiksaanwijzing
Instruções de Utilização
Οδηγίες χρήσης
Brugervejledning
Bruksanvisning
Käyttöohjeet
Bruksanvisning
Instrukcja użycia
Használati útmutató
Návod k použití
Kullanım Talimatları
使用説明
사용 설명서
Инструкции по Применению

Single Use

Nonabsorbable

Rx only

STERILE EO

Do not resterilize

BARD

DAVOL INC.

BARD® Mesh Pre-shaped

Monofilament Polypropylene for the Reconstruction of Tissue Defects

INSTRUCTIONS FOR USE

Product Description:

BARD® Pre-shaped mesh is a preshaped mesh constructed of knitted polypropylene monofilaments. The knit construction allows the mesh to be stretched in both directions, to accommodate and reinforce tissue defects. BARD® Pre-shaped mesh is a nonabsorbable, sterile prosthesis.

Indications:

BARD® Pre-shaped mesh is indicated for the repair of inguinal hernia defects.

Contraindications:

1. Do not use BARD® Pre-shaped mesh in infants or children, whereby future growth will be compromised by use of such mesh material.
2. Literature reports there may be a possibility for adhesion formation when BARD® mesh is placed in direct contact with the bowel or viscera.

Warnings:

1. This device must be sterile before use. Please inspect the packaging to be sure it is intact and undamaged.
2. This device is for single use only. Do not resterilize or reuse any portion of BARD® Mesh Pre-shaped.
3. Careful attention to BARD® Pre-shaped mesh handling, fixation, and suture technique is most important in the presence of known or suspected wound contamination or infection.
4. The use of any permanent mesh or patch in a contaminated or infected wound can lead to fistula formation and/or extrusion of the prosthesis.
5. If an infection develops, treat the infection aggressively. Consideration should be given regarding the need to remove the mesh. An unresolved infection may require removal of the device.
6. To prevent recurrences when repairing inguinal hernias, the mesh should be large enough to extend beyond the pubic tubercle and should fit securely around the spermatic cord at the internal ring. Many surgeons cut a keyhole in the mesh to allow for easier placement around the cord.¹
7. Reuse, reprocessing, reesterilization or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device and may lead to device failure which may result in injury to the patient.

Reuse, reprocessing, reesterilization or repackaging may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient or end user.

8. If unused prosthesis has been in contact with instruments or supplies used on a patient or contaminated with body fluids, discard with care to prevent risk of transmission of infection.

Precautions:

1. Please read all instructions prior to use.
2. Only physicians qualified in the appropriate surgical techniques should use this prosthesis.
3. Intact BARD® Pre-shaped mesh exhibits high burst and tensile strength. However, when custom tailoring, in special circumstances where excessive force is placed on the mesh, the following guidelines may be helpful:
 - When cutting a notch in the mesh, a V-shape with a radiused point will withstand more force than a V-shape which comes to a sharp point.
4. DAVOL™ permanent or absorbable fixation devices or nonabsorbable monofilament sutures are recommended to properly secure the prosthesis. If absorbable fixation devices are used, they must be indicated for use in hernia repair.
5. Care should be taken to ensure that the mesh is adequately fixated to the uncompromised tissue of the inguinal floor. If necessary, additional fasteners and/or sutures should be used.

Adverse Reactions:

Possible complications include seroma, adhesions, hematoma, inflammation, extrusion, fistula formation and recurrence of the hernia or soft tissue defect.

Fixation:

Care should be taken to ensure that the mesh is adequately fixated to the uncompromised tissue of the inguinal floor. If necessary, additional fasteners and/or sutures should be used.

DAVOL™ permanent or absorbable fixation devices or nonabsorbable monofilament sutures are recommended to properly secure the prosthesis. If absorbable fixation devices are used, they must be indicated for use in hernia repair. The minimum distance for fixation points should be 4 mm from the mesh edge.

Traceability:

A traceability label is attached to every prosthesis package which identifies the type, size, and lot number of the prosthesis. This label should be affixed to the patient's permanent medical record to clearly identify the device which was implanted.

Copyright © 2005-2008, 2011, 2013 C. R. Bard, Inc.
All Rights Reserved.

Bard and Davol are trademarks and/or registered trademarks of
C. R. Bard, Inc. or an affiliate.

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1. Lichtenstein, I.L., Hernia Repair Without Disability 2nd Edition St. Louis, 1986, Ishiyaku EuroAmerica, Inc.
2. Martin, R.E., et al., "Polypropylene Mesh in 450 Hernia Repairs: Evaluation of Wound Infections", *Contemporary Surgery*, 1982 20:46-48.
3. Barnes, J.P., "Inguinal Hernia Repair with Routine Use of Marlex® Mesh", *Surgery, Gynecology and Obstetrics*, 1987 165:33-37.
4. Lichtenstein, I.L., et al., "Cause and Prevention of Post Herniorrhaphy Neuralgia: A Proposed Protocol for Treatment", *The American Journal of Surgery*, 1988 155:786-790.
5. Lichtenstein, I.L., "Herniorrhaphy", *The American Journal of Surgery*, 1987 153:553-559.
6. Martin, R.E., Max, C.C., "The Use of Marlex® Mesh in Primary Hernia Repairs", *Surgical Rounds*, 1983 6:52-62.

Manufacturer:

Davol Inc.
Subsidiary of C. R. Bard, Inc.
100 Crossings Boulevard
Warwick, RI 02886 USA
1-401-825-8300 • 1-800-556-6275

Bard Limited
Crawley, UK
RH11, 9BP

Medical Services & Support
Clinical Information Line
1-800-562-0027

PK3797915
142R

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Copyright © 2005-2008, 2011, 2013 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2005–2008, 2011, 2013 C. R. Bard, Inc.
Copyright © 2005-2008, 2011, 2013 C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2005-2008, 2011, 2013 C. R. Bard, Inc. Todos os direitos reservados.
Copyright © 2005-2008, 2011, 2013 C. R. Bard, Inc. Alle rettigheder forbeholdes.
Copyright © 2005-2008, 2011, 2013 C. R. Bard, Inc. Kaikki oikeudet pidätetään.