PROMAX™ Series

BARD® V-MAX™ 24 Roller Bar

Instructions for Use

Contents:

- 2 - Compatibility Table
- 3 - Symbols
- 6 - English
- 8 - French/Français
- 10 - Spanish/Español

BARD

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* = 0.014" wire size.

Manufacturer
Fabricant
Fabricante

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R_x Only

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practices and applicable local, state and federal laws and regulations.

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ProMAX™ Series**BARD® V-MAX™24 Roller Bar****Information For Use****Device Description:**

The BARD® V-MAX™24 Roller Bar device is a monopolar electrode designed to deliver radio frequency energy that is supplied by an electrosurgical generator cleared for medical use.

Indications For Use:

The BARD® V-MAX™24 Roller Bar device is indicated for ablation and coagulation of soft tissue including procedures for the treatment of benign prostatic hypertrophy (BPH). It is intended for use with compatible resectoscopes.

Potential Complications:

Complications could include local and/or systemic infection, thermal damage to surrounding structures, local hematoma, dissection and perforation, distal tip separation, and patient discomfort during and/or after energy application.

Precautions:

1. All medical staff should carefully review product labeling and instruction sheets before using the BARD® V-MAX™24 Roller Bar device. Inappropriate use of the instrument could adversely affect the procedure or cause injury to the patient.
2. Refer to the applicable operating and maintenance manuals for the resectoscope and electrosurgical generator being used.
3. The BARD® V-MAX™24 Roller Bar device should be used only by a physician who is familiar with the use of electrosurgical instruments, devices and power generators. Consult the medical literature regarding techniques, typical power settings, complications, and hazards prior to any endoscopic procedure.
4. During application of electrical energy, movement of the roller is required to achieve the desired tissue effect.
5. The time and energy required for treating specific tissue may differ when using the BARD® V-MAX™24 Roller Bar device compared to other electrosurgical devices.
6. The monopolar electrosurgical generator should be set to the "Cut" or "Pure Cut" mode for maximum tissue removal effect. Do not use the "Blend" mode and be cautious of arcing on high coagulation power settings.
7. Starting at power settings consistent with standard resectoscope cutting loop procedures, gradually increase power until the desired tissue effect is achieved. If there is little or no tissue effect, check the generator, power and grounding cables and the electrodes. Various electrosurgical generators have different power delivery capabilities. Depending on the impedance in the monopolar circuit, the power delivered may be lower than the power setting indicates (see owner/operator manual for the generator being used).
8. Constant irrigation is required throughout the procedure, and the distal tip of the device should be submerged and kept in view at all times. Use sterile non-conductive irrigation solution only.
9. Immediately discontinue use if breaks or fractures appear in the BARD® V-MAX™24 Roller Bar device. Breaks or fractures may allow undirected emission of electrical energy, rendering the device useless and potentially causing harm to surrounding tissues.
10. **Do not bend or manipulate the device.**
11. Care should be taken to avoid severe impacts, side stresses or bends at sharp angles.
12. When endoscopic devices are used together, ensure that any isolation or ground is not violated.
13. If a standard resectoscope cutting loop is used after using the BARD® V-MAX™24 Roller Bar device, the electrosurgical power setting must be adjusted to the appropriate level.

How Supplied and Stored:

- The BARD® V-MAX™24 Roller Bar device is supplied sterile by EtO and is designed for single use.
- Store in a cool, dry place.
- Sterile unless package is opened or damaged.

Directions For Use:

1. Remove the BARD® V-Max™24 Roller Bar device from the package and examine it for damage. Do not use if there are visible signs of damage or insulation is not intact.
2. Assemble the working element and resectoscope according to the instructions provided by the resectoscope manufacturer.
3. Insert the BARD® V-Max™24 Roller Bar device into the working channel of the resectoscope in the same manner as for an existing cutting loop electrode and consistent with the instructions of the resectoscope.
4. Ensure that the electrode is securely placed and locked into position by pulling carefully on the electrode stabilizer sleeve (if appropriate).
5. Attach resectoscope to the electrosurgical generator according to manufacturer's recommendations.
6. Insert the combined assembly into the resectoscope sheath and position at the point where initial application of electrical energy will be delivered.
7. Unless recommended by other sources, the electrosurgical generator should be set on the "Pure Cut" setting at power levels consistent with standard resectoscope roller bar procedures.
8. The power level can be adjusted and varied slightly for different patients and also when tissue is being removed on subsequent sweeps.
9. The level of effect is governed by the amount of time the electrode is in contact with an area and the power delivered.
10. The slower the electrode is moved across any particular area and the higher the power applied, the deeper the level of effect.
11. Maintain constant irrigation through the resectoscope using sterile non-conductive irrigant. Keep the roller in the field of view and submerged in irrigant at all times.
12. If tissue adheres to the roller bar, cleaning can be accomplished by switching to a coagulation setting and engaging current without tissue contact and with roller fully submerged in irrigation solution. Cleaning can also be accomplished by removing the device and wiping it with a soft, sterile lint-free cloth.

Bard Medical Division warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Medical Division sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD MEDICAL DIVISION BE LIABLE TO YOU FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

Customer Service:

In the event the device must be returned for any reason, contact C. R. Bard, Inc. Customer Service at (800) 526-4455 if inside the U.S., to receive a return authorization number prior to return shipment. If outside the U.S., contact your local Customer Service representative. Return in original package and shipping box. For more information about the BARD® V-Max™24 Roller Bar device or other BARD® products, contact Customer Service or your local Customer Service representative.

This is a single use device. Do not resterilize any portion of this device. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.

Manufactured in U.S.A. • Packaged in Mexico.

Bard, ProMax and V-Max are trademarks and/or registered trademarks of C. R. Bard, Inc.
©2015 C. R. Bard, Inc. All rights reserved.

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Bard, ProMax y V-Max son marcas comerciales y/o comerciales registradas de C. R. Bard, Inc.
©2015 C. R. Bard, Inc. Reservados todos los derechos.