

PRESCRIPTIVE INFORMATION

Indications for Use: This device is intended for use in the endoscopic removal of renal and ureteral stones.

Warnings

- Some objects may be too large to be removed endoscopically using a retrieval device. The use of fluoroscopy and/or x-ray to determine the size of the object is recommended; do not use the stone basket if the object is too large to be removed endoscopically.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Caution

- Objects that are too large to be recovered through the sheath or through the scope channel will require the scope and the basket to be removed simultaneously from the urinary tract.
- If resistance is encountered during advancement or withdrawal of the device, stop and determine the source of the resistance, as continued resistance may damage the device and could result in patient injury. Take action to alleviate the resistance. Where necessary, use of the lithotripsy may be required to reduce the stone burden within the basket, providing that no direct contact is made with the stone basket.
- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Precautions

- Do not allow the device to come in contact with any electrified instrument.
- Do not allow the device to be directly fired upon by any lithotripsy device. To do so could damage the device and could result in patient injury.

Potential Complications that may result from the use of a basket in an endoscopic urological procedure include:

Perforation	Evulsion	Edema	Entrapment
Basket Inversion	Hemorrhage	Inability to disengage from irretrievable object	

Notes

- Only physicians trained in stone management should perform this procedure. A variety of techniques in the use of this instrument may be employed; however the physician should use the technique with which they are most familiar with.
- Single patient use, do not reuse.
- Sterile unless package is opened or damaged.
- Do not re-sterilize.