TRUE[®] FLOW

Valvuloplasty Perfusion Catheter

ENGLISH

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Device Description

The TRUE[®] FLOW Valvuloplasty Perfusion Catheter is an over-the-wire co-axial catheter with a balloon fixed at the tip. The balloon enables hemodynamic flow through its central orifice. The catheter is 110cm long and has two lumens: one lumen is used to inflate and deflate the balloon and the other permits the use of a guidewire to position the catheter. The balloon inflation luer-lock hub (angled) connects to a syringe inflation device to deliver radiopaque contrast media for inflation. The guidewire luer-lock hub (straight) connects to the guidewire lumen. The balloon is non-compliant and is designed to reach a known diameter and length when inflated within the specified pressure range. Two radiopaque markers are located on the guidewire lumen. These bands are positioned at the proximal and distal balloon shoulders. These markers are provided for fluoroscopic positioning of the device across the aortic valve. Balloon catheter dimensions, nominal pressure, maximum inflation pressure, recommended introducer size, and maximum guidewire size are indicated on the package label.

Indications for Use

The TRUE[®] FLOW Valvuloplasty Perfusion Catheter is indicated for balloon aortic valvuloplasty.

Contraindications

None known

Warnings

1. Do not use in patients with annular dimensions <18mm.
2. Contents supplied STERILE using ethylene oxide (EO). Non-pyrogenic. Do not use if sterile barrier is opened or damaged. Single patient use only. Do not reuse, reprocess or re-sterilize.
3. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminate amount of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
4. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminate degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
5. Catheter balloon inflation diameter must be carefully considered in selecting a particular size for any patient. It is critical to perform a clinical diagnostic determination of valve anatomical dimensions prior to use; imaging modalities such as transthoracic echocardiogram (TTE), computerized tomography (CT), angiography, and/or transesophageal echocardiogram (TEE) should be considered. The inflated balloon diameter should not be significantly greater than valvular diameter.
6. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip breakage or balloon separation, or cause injury to the patient (such as vessel perforation).
7. If flow through catheter becomes restricted, do not attempt to clear catheter lumen by infusion. Doing so may cause catheter to rupture, resulting in vessel trauma. Remove and replace catheter.
8. Do not exceed the RBP recommended for this device. Balloon rupture may occur if the RBP rating is exceeded. To prevent over-pressurization, use of a pressure monitoring device is recommended.
9. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state, and federal laws and regulations.
10. If using device to support Transcatheter Aortic Valve Implantation (TAVI), consult TAVI system's Instructions for Use for any additional procedural instructions related to selection and use of valvuloplasty balloon.
11. To reduce thrombosis, this device should not be used without appropriate anticoagulation. It is recommended to maintain an ACT of ≥ 200 seconds during use of this device.

Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. The catheter should only be used by physicians trained in the performance of percutaneous transluminal valvuloplasty.
3. The minimal acceptable French size is printed on the package label. Do not attempt to pass the catheter through a smaller size sheath introducer than indicated on the label.
4. Use the recommended balloon inflation medium of 1/3 to 2/3 contrast to saline ratio. Never use air or other gaseous medium to inflate the balloon.
5. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath, determine if contrast is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the sheath and then completely evacuate the contrast before proceeding to withdraw the balloon.
6. If resistance is still felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and guidewire/introducer sheath as a single unit.
7. In the very unlikely event of balloon burst or rupture, balloon could be more difficult to remove through the sheath and could require introducer sheath removal.
8. Do not torque, excessively bend catheter or continue to use if the shaft has been bent or kinked.
9. Do not re-insert the catheter into the body once it has been removed from the sheath, as withdrawing balloon through introducer sheath may damage balloon.
10. Do not remove guidewire from catheter during procedure.
11. Dilation procedures should be conducted under high-quality fluoroscopic guidance.
12. Careful attention must be paid to the maintenance of tight catheter connections. Aspirate before proceeding to avoid air introduction into the system.
13. If inflating balloon in patient to facilitate re-folding, ensure balloon is positioned so that it can be inflated safely.

Potential Adverse Reactions

The complications which may result from a percutaneous transluminal valvuloplasty procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Cardiovascular injury
- Conduction system injury
- Death
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stroke
- Thrombosis
- Valvular tearing or trauma
- Vessel dissection, perforation, rupture, or spasm

Directions for Use

Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources.

Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date.

Do not use if packaging is damaged or opened.

Equipment for Use

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (50 ml or larger)
- Appropriate introducer sheath and dilator set
- .035" guidewire
- 3-way high pressure stopcock

Dilatation Catheter Preparation

1. Remove catheter from package. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
2. Keep balloon guard on balloon until removing it immediately prior to inserting catheter into introducer.
3. Mix contrast / saline with a ratio of 1/3 contrast, 2/3 saline by volume.
4. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 50ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium. Do not use air or any gaseous medium to inflate the balloon.
5. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
6. Connect the syringe to the stopcock.
7. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
8. Repeat step #7 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
9. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.

Use of the TRUE® FLOW Valvuloplasty Perfusion Catheter

1. Prepare patient and place introducer and guidewire into access site using standard percutaneous catheterization techniques.
2. Remove balloon guard distally, ensuring guard is no longer on catheter.
3. Backload the distal tip of the TRUE® FLOW Valvuloplasty Perfusion Catheter over the prepositioned guidewire and advance the tip to the introduction site.
4. Using fluoroscopic guidance, advance the catheter through the introducer sheath and over the wire to the site of inflation. Discontinue product use if marker bands are not visible under fluoroscopy.

5. Position the balloon within the relevant area of the aortic valve to be dilated, ensure the guidewire is in place, and while ensuring the balloon is held in a static position, inflate the balloon to a pressure not greater than RBP.

Note: Blood will flow through the balloon.


6. Apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy. If balloon does not appear to be fully deflated, ensure product is outside aortic valve in safe position for secondary inflation. Inflate slightly and deflate again to ensure complete deflation is achieved.
7. While maintaining negative pressure and the position of the guidewire, withdraw deflated catheter over the guidewire and through the introducer sheath. Use of a gentle clockwise twisting motion may be used to help facilitate catheter removal through the introducer sheath.
8. If unusual resistance is met when attempting to withdraw balloon, position balloon in anatomical position in which it can be inflated safely. Inflate balloon and then deflate it. Balloon re-folding can be observed under fluoroscopy. Use of recommended contrast concentration will facilitate fluoroscopic visualization of balloon re-folding.

Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet.

In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

Technical Specifications

Simulated Hemodynamic Testing Results*

Measure (Averages)	Result
Cardiac Output	1.0 L/min
Mean Gradient	74.7 mmHg
Aortic Regurgitation	74.6%
Aortic Systolic Pressure	50.2 mmHg
Aortic Diastolic Pressure	7.7 mmHg
Aortic Mean Pressure	26.0 mmHg

* Bench testing conducted using a Vivitro Pulse Duplicator System on an 18 mm inflated balloon, heart rate of 70 bpm, systolic waveform of 35%, cardiac output of 3.5 L/min and pressure of 120/80 mmHg. N=30. Bench data results may not be indicative of actual clinical performance. Different tests may yield different results.

Catalogue No.		TF0183511		TF0203511		TF0223512		TF0243514		TF0263516	
Balloon Description		18 mm x 3.5 cm		20 mm x 3.5 cm		22 mm x 3.5 cm		24 mm x 3.5 cm		26 mm x 3.5 cm	
Actual Diameter (mm) Over Indicated Pressure Range (3 atm - 6 atm)	Pressure	3 atm	6 atm	3 atm	6 atm	3 atm	6 atm	3 atm	6 atm	3 atm	6 atm
	Dim A	18.4	18.4	20.3	20.4	22.3	22.4	24.6	24.7	26.4	26.6
	Dim B	17.5	17.7	19.4	19.5	21.3	21.4	23.5	23.7	25.3	25.5

Note: The indicated dimensions are based on bench top testing of an unconstrained balloon at operating pressure.

Clinical Experience

Design

The primary objective of the TRUE-FLOW study was to observe the performance and safety of the TRUE® FLOW Valvuloplasty Perfusion Catheter during dilatation of the aortic valve in preparation for Transcatheter Aortic Valve Implantation (TAVI). Twenty-five subjects were enrolled, of which all were successfully treated with the study device. One subject withdrew consent following all study procedures and hospital discharge, and is thus not included in the Analysis Population (n = 24).

Subjects eligible to be enrolled in this study were adult male and non-pregnant female patients scheduled to undergo TAVI for the treatment of aortic stenosis who had an annulus diameter that could be treated with the available sizes of the study device, in accordance with the IFU.

Demographics, Medical History & Risk Calculations

The mean age for patients participating in the study was 79.8 years old, and 70.8% of patients were male. The mean Society of Thoracic Surgeons (STS) score was 7.8%. COPD was present in 79.2% of patients, atrial fibrillation in 33.3%, coronary artery disease in 87.5% and peripheral artery disease in 58.3%. 45.8% of patients were diabetic while 4.2% of patients were considered to be in renal failure. All patients were considered to be in congestive heart failure with a mean (SD) ejection fraction of 53.1% (13.87%).

With regards to prior interventions, no patients had received prior surgical valve treatment or balloon aortic valvuloplasty (BAV). One patient (4.2%) had prior coronary artery bypass graft surgery while six (25%) had undergone a percutaneous coronary intervention. Three patients (12.5%) had pre-existing pacemakers/defibrillators.

The mean (SD) aortic valve area was 485.6mm² (88.73mm²) with an average annulus diameter of 25.4mm (2.18mm). Mean (SD) annulus calcification grade was 267.5mm³ (300.84mm³) using a modified Agatston methodology for contrast CT scans and a detection threshold of 850 Hounsfield Units.

Aortic Valve Characteristics (n=24)	
Characteristic	Mean ± SD (Min, Max)
Aortic Valve Area (mm ²)	485.6 ± 88.73 (307, 614)
Perimeter Derived Aortic Valve Area (mm ²)	516.8 ± 92.79 (327, 642)
Annulus Diameter (mm)	25.4 ± 2.18 (20, 29)
Annulus Calcification Grade (mm ³)*	267.5 ± 300.84 (1,1082)

*Annulus Calcification Grade: Determined using modified Agatston methodology for contrast scans and a detection threshold of 850 Hounsfield units.


The mean (SD) STS score was 7.8% (5.87%) for the patients enrolled in the trial. Thirteen (54.2%) patients were in Class III with eight (33.3%) patients in Class II and 3(12.5%) of patients in Class IV.

Risk Calculations (n=24)	
Characteristic	Mean ± SD (Min, Max)
STS Score (%)	7.8 ± 5.87 (2, 25)
EURO Score I (%)	16.3 ± 11.03 (5, 53)
EURO Score II (%)	8.3 ± 7.50 (2, 25)
NYHA Classification	n (%)
Class II	8 (33.3)
Class III	13 (54.2)
Class IV	3 (12.5)

Results


In 21 (91.3%) cases the TRUE® FLOW catheter successfully dilated the aortic annulus without clinically significant movement based on a qualitative judgement by the operator. For one subject these data were not recorded. On average (SD) the device was delivered within 4.2mm (5.0mm) of the target position and with an average (SD) of 2.0mm (4.5mm) of movement while the device was fully inflated. No ventricular pacing was used in 22 (91.7%) of the study cases. Two patients, both with preexisting permanent pacemakers, were moderately paced at 120 BPM and 160 BPM, respectively. This study was conducted in a patient population with cardiovascular risk factors known to be associated with increased occurrence of complications following rapid pacing (RP) during BAV.

The mean (SD) total time of balloon inflation was 12.5 seconds (3.74 seconds). The balloon remained fully inflated for a mean (SD) of 3.2 seconds (1.06 seconds).

No device related serious adverse events were reported during the study. One patient developed a case of procedure related pericardial effusion before the introduction of the study device which resolved with appropriate treatment. No other serious adverse events or unanticipated device effects occurred during the study.

Additionally, intraventricular pressure measurements using the transeptal approach were taken before, during, and after the index procedure and were used to assess the ability of the inflated TRUE® FLOW catheter to remain stationary with clinically acceptable intraventricular pressure under rapid pacing, reduced pacing, or without pacing. For all 24 procedures, mean (SD) left ventricular end diastolic pressure (LVEDP) during the procedures [13.4 (11.33) mmHg] was lower than that prior to the procedures [17.8 (10.88) mmHg].

The results of the TRUE-FLOW prospective, observational study indicated that the design of the TRUE® FLOW Valvuloplasty Perfusion Catheter allowed adequate blood flow through the device's central orifice to perform complete dilation of a stenotic aortic valve prior to TAVI without the need for rapid pacing.

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Copyright ©2018 C. R. Bard, Inc. All rights reserved. Printed in the USA.
Copyright ©2018 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright ©2018 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in USA.
Copyright ©2018 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright ©2018, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
Copyright ©2018 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.