DISTAFLLO[®]

Bypass Grafts

Instructions for Use

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INSTRUCTIONS FOR USE

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Device Description, Indications, Contraindications, Warnings, Precautions and Adverse Reactions

Device Description

DISTAFLO® bypass grafts are constructed of expanded polytetrafluoroethylene (ePTFE) with carbon impregnated into the inner portions of the graft wall. This graft is also modified with a proprietary cuff at the distal end.

DISTAFLO® bypass grafts have a removable external spiral bead support over the full length of the graft and can be used where resistance to compression or kinking is desired.

DISTAFLO® bypass grafts with standard and Small Cuff have a one (1) mm internal diameter reduction over the distal 25 cm of the graft. The DISTAFLO® Mini-Cuff graft has a two (2) mm internal diameter reduction over the distal 25 cm of the graft. This reduction is intended to increase the velocity of flow as it approaches the cuff area.

Indications for Use

DISTAFLO® bypass grafts are intended for bypass or reconstruction of peripheral arterial blood vessels.

Contraindications

None known.

Warnings

1. All DISTAFLO® bypass grafts are supplied sterile and non-pyrogenic unless the package is opened or damaged. DISTAFLO® bypass grafts are sterilized by ethylene oxide.
2. Do not use after expiration date printed on the label.
3. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
4. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
5. Anastomotic or graft disruption has been associated with Axillofemoral, Femoral Femoral, or Axillobifemoral bypass procedures if implanted improperly. Refer to *Specific Operative Procedures (Extra-Anatomic Bypass Procedures)* for further instructions.
6. For Extra Anatomic procedures (e.g. Axillofemoral, Femoral Femoral, or Axillobifemoral bypass), the patient should be cautioned that sudden, extreme or strenuous movements should be totally avoided for a period of at least six to eight weeks to allow for proper stabilization of the graft. Routine activities such as raising the arms above the shoulders, reaching out in front, extended reaching, throwing, pulling, striding or twisting should be avoided.
7. DISTAFLO® bypass grafts do not stretch (are non-elastic) in the longitudinal direction. The correct graft length for each procedure must be determined by considering the patient's body weight, posture, and the range of motions across the anatomical area of graft implantation. Failure to cut the graft to an appropriate length may result in anastomotic or graft disruption, leading to excessive bleeding, and loss of limb or limb function, and/or death.
8. Aggressive and/or excessive graft manipulation when tunneling, or placement within a too tight or too small tunnel, may lead to separation of the spiral beading and/or graft breakage. The distal anastomosis should be made after tunneling or suture disruption can occur. DO NOT pass the cuff portion (distal end) of the DISTAFLO® bypass graft through a tunneler sheath or the tissue tunnel, as this could lead to separation of the spiral beading and/or graft breakage.
9. When embolectomy or balloon angioplasty catheters are used within the lumen of the graft, the inflated balloon size must match the inner diameter of the graft. Over-inflation of the balloon or use of an inappropriately sized balloon may dilate or damage the graft.
10. Avoid repeated or excessive clamping at the same location on the graft. If clamping is necessary, use only atraumatic or appropriate vascular smooth jawed clamps to avoid damage to the graft wall. Do not clamp the cuffed portion of the graft.

11. Exposure to solutions (e.g., alcohol, oil, aqueous solutions, etc.) may result in loss of the graft's hydrophobic properties. Loss of the hydrophobic barrier may result in graft wall leakage. Preclotting of this graft is unnecessary.
12. Avoid excessive graft manipulation after exposure to blood or body fluids. Do not forcibly inject any solution through the lumen of the graft, or fill the graft with fluid prior to pulling it through the tunnel as loss of the graft's hydrophobic properties may occur. Loss of the hydrophobic barrier may result in graft wall leakage.
13. Do NOT expose DISTAFLO® bypass grafts to temperatures greater than 500°F (260°C). PTFE decomposes at elevated temperatures, producing highly toxic decomposition products.¹
14. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.
15. During tunneling, create a tunnel that closely approximates the outer diameter of the graft. A tunnel that is too loose may result in delayed healing and may also lead to perigraft seroma formation.

Precautions

1. Only physicians qualified in vascular surgery techniques should use this prosthesis. The healthcare provider is responsible for all appropriate postoperative care instructions to the patient.
2. The healthcare provider must observe aseptic technique during implantation and postoperatively.
3. When removing the external spiral support (beading) of the DISTAFLO® graft, the beading must be removed slowly and at a 90° angle to the graft. Rapid unwinding and/or removal at less than a 90° angle may result in graft damage. Do not use surgical blades or sharp, pointed instruments to remove the beading as this may damage the graft wall. If damage occurs, that segment of the graft should not be used. Refer to "Anastomotic Preparation", for further instructions. (Reference Figure 1.)

Figure 1

4. DISTAFLO® grafts have been developed for and are especially suitable for below the knee and infrapopliteal bypass and are not recommended for Extra-Anatomic Bypass Applications.
5. When suturing, avoid excessive tension on the suture line, inappropriate suture spacing and bites, and gaps between the graft and host vessel. Failure to follow correct suturing techniques may result in suture hole elongation, suture pull-out, anastomotic bleeding and/or disruption. Refer to "Suturing" for further instructions.
6. To minimize fluid collection around the graft in Extra-Anatomic bypass procedures or in peripheral reconstruction procedures, the lymphatics should be carefully ligated and sealed, especially in the groin area.
7. Consider intraoperative and postoperative patient anticoagulation therapy for each patient as appropriate.

Adverse Reactions

Potential complications which may occur with any surgical procedure involving a vascular prosthesis include, but are not limited to: Disruption or tearing of the suture line, graft, and/or host vessel; suture hole bleeding; graft redundancy; thrombosis; embolic events; occlusion or stenosis; ultrafiltration; seroma formation; swelling of the implanted limb; formation of hematomas or pseudoaneurysm; infection; skin erosion; aneurysm/dilation; blood leakage; and hemorrhage.

DIRECTIONS FOR USE

Equipment Required

Tunneler, double armed sutures, atraumatic clamps, scalpels and/or scissors.

Opening the Package

Hold the outer tray in one hand. Peel back the lid. Remove the inner tray. Peel back the inner tray lid slowly and carefully remove the graft using sterile atraumatic instruments or gloves. Protect the graft against damage from sharp or heavy instruments.

General Operative Techniques

Sizing

DO NOT MODIFY THE DISTAFLO® CUFF.

DISTAFLO® grafts are available in standard, Small and Mini-Cuff configurations. The Small Cuff configuration should be used in recipient vessels with diameters of 4mm or smaller, only, and the Mini-Cuff configuration should be used in recipient vessels with diameters of 3 mm or smaller only.

Tunneling Recommendations

Note: Prior to utilizing a sheath tunneler, verify that the DISTAFLO® graft O.D. dimensions fit the sheath I.D. dimensions.

Use of a sheath tunneler is recommended for placing the DISTAFLO® bypass graft. Use of a sheath tunneler will minimize graft handling, protect the external reinforcement (Flex beading) and help maintain graft integrity.

Always follow the instructions for use for the specific tunneler utilized to place the DISTAFLO® bypass graft.

Create a tunnel that closely approximates the outer diameter of the graft. A tunnel that is too loose may result in delayed healing and may also lead to perigraft seroma formation. **See Warnings #8 and #15.**

Anastomotic Preparation of DISTAFLO® Grafts with Flex Beading

Use sterile technique to remove the Flex beading. Hold the graft firmly with one gloved hand. With the other gloved hand, gently grasp and lift the end of the beading with gloved fingers. Slowly unwind the beading by peeling at a 90° angle to the graft (**reference Figure 1 and Precaution #3**). Remove and trim off enough beading to cut the graft to appropriate length at the non-cuffed end.

Suturing

Size the graft appropriately to minimize excessive tension at the suture line. Use a tapered, non-cutting needle with a nonabsorbable monofilament suture approximately the same size as the needle. Take 2 mm suture bites in the graft following the curve of the needle and gently pull the suture at a 90° angle. Proper sizing of the graft length prior to implant will minimize suture hole elongation caused by excessive tension. **Reference Warning #7 and Precaution #5.**

Suggested Suturing Technique


Perform an arteriotomy five to eight millimeters shorter than the length of the cuff.

Using double armed suture, tack the heel of the cuff, aligning the toe in the direction of distal blood flow.

Suture continuously from the heel to just beyond the midpoint of the cuff on one side of the arteriotomy. Using the other arm of the heel suture, continue suturing to just beyond the midpoint of the graft on the other side of the arteriotomy.


Extend the arteriotomy to the exact cuff length (**reference Figure 2**).

Figure 2

Complete the "cuffed" anastomosis by suturing with one of the suture arms, around the toe, to the other side of the cuff. Always tie off on the side of the cuff (**see Figure 3**).

Figure 3

Perform the proximal arterial anastomosis using standard surgical techniques. In order to retain the cuffed configuration, only the proximal end of the DISTAFLO® graft should be sized for length.

Thrombectomy

Techniques for declotting DISTAFLO® bypass grafts include but are not limited to the use of balloon catheters. **See Warning #9.**

Longitudinal Incision:

Place stay sutures before introducing the embolectomy catheter. Cut through the spiral beading and base graft. The spiral beading will realign itself after closure. Place a longitudinal incision in the graft that is long enough to accommodate the extraction of a fully dilated thrombectomy catheter balloon. A patch may be considered as an aid to graft closure.

Transverse Incision:

No stay sutures are necessary. A horizontal mattress suture is recommended for graft closure.

During the early postoperative period, the natural progression of wound healing renders the graft translucent in appearance. In this state, a longitudinal incision with stay sutures is recommended. If a transverse incision is performed, a horizontal mattress suture technique and PTFE pledgets may aid in closure.

Angiography

Should angiography be performed at the time of the procedure, the artery proximal to the graft should be used for injection, if possible.

Specific Operative Procedures

Extra-Anatomic Bypass Procedures (e.g., Axillofemoral, Femoral Femoral, and Axillofemoral)

For Extra-Anatomic bypass procedures, careful attention must be given to the following techniques. Failure to follow these technical considerations may result in suture hole elongation, mechanical disruption, or tearing of the graft, suture line or host vessel, thrombosis, extreme blood loss, loss of limb function, loss of limb, or death. (**Reference Precaution #4 and #5 and Warning #6 and #7.**)

- To avoid extreme stress on the anastomosis and the graft, include the patient's weight and range of limb motion when determining graft length, tunnel length and location.
- To determine the correct graft length, drape the patient to allow full movement of the arm, shoulder girdle or legs.
- Avoid protracted hyperabduction of the arm, during the surgical procedure. Prolonged hyperabduction may lead to brachial plexus injury.
- Allow sufficient graft length to avoid stressing of axillary or femoral anastomosis throughout the full range of movement of the arm, shoulder girdle, or legs. The graft should be placed under both the pectoralis major and pectoralis minor. (**Reference Figure 4.**)

Figure 4

- Cutting the graft slightly longer than necessary has been reported by some surgeons to further reduce the risk of stressing the graft or the anastomosis.
- Correctly bevel the axillary anastomosis. Stress on the graft is minimized when the graft is placed perpendicular (0°) to the axillary artery. Therefore, the anastomotic angle should be as small as possible and should not exceed 25° relative to the cut edge of the graft.
- Place the graft anastomosis close to the rib cage on the first portion of the axillary artery.² Do not place the anastomosis on the third portion of the axillary artery.
- Notify the patient that sudden, extreme or strenuous movements of the arm, shoulder or leg should be totally avoided for a period of at least six to eight weeks to allow for proper stabilization of the graft. Routine activities such as raising the arms above the shoulders, reaching out in front, extended reaching, throwing, pulling, striding, or twisting should be avoided. **See Warning #5, #6 and #7, and Precaution #5.**

References

- Guide to Safe Handling of Fluoropolymer Resins, Fourth Edition, The Fluoropolymer Division of the Society of Plastics Industry, Inc.
- Victor M. Bernhard, M.D. and Jonathan B. Towne, M.D., Editors, *Complications in Vascular Surgery*, Second Edition, Grune and Stratton, Inc. (Harcourt Brace Jovanovich; Publishers), Orlando, 1985, 56.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Copyright © 2016 C. R. Bard, Inc. All rights reserved. Printed in the U.S.A.**
Copyright © 2016 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright © 2016 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright © 2016 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright © 2016, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
Copyright © 2016 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.