# ULTRAScore™

## Focused Force PTA Balloon

ENGLISH

### INSTRUCTIONS FOR USE

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

#### Device Description

The ULTRAScore™ Focused Force PTA Balloon consists of a flexible, over-the-wire (OTW) catheter shaft with a semi-compliant balloon fixed at the distal end. For all balloon lengths, radiopaque markers delineate the working length of the balloon and aid in balloon placement. For balloon lengths of 100 mm and greater, two radiopaque markers are positioned on the distal portion of the balloon and one radiopaque marker is positioned on the proximal portion of the balloon to differentiate between the distal and proximal ends of the balloon. The catheter includes an atraumatic tip to facilitate advancement of the catheter to and through the stenosis. Two scoring wires, oriented 180° apart, provide focused force upon dilatation. The ULTRAScore™ Focused Force PTA Balloon is compatible with .014" or .035" guidewires, as denoted by the product labeling. The distal portion of the .014" guidewire compatible catheters is hydrophilically coated. The proximal portion of the catheter includes a female luer lock hub connected to the catheter with a guidewire lumen and an inflation lumen. Packaged with every product is a protective sheath that is positioned over the balloon and must be removed prior to use. A stylet is placed into the tip of the catheter. These products are not made with natural rubber latex.

Figure 1

The GEOALIGN® Marking System is a non-radiopaque ruler on the catheter shaft measured from the distal tip. The GEOALIGN® markings are designated on the catheter shaft by 1cm increment bands with an accuracy within ±1mm. The distance from the distal catheter tip is labeled in 10cm increments. Thicker bands denote the midway point (5cm) between the labeled distances. The GEOALIGN® Marking System is designed to be used as a tool to externally measure the intravascular advancement and/or retraction of the catheter. This can provide an intravascular reference regarding the location of the distal tip of the catheter or an approximate intravascular length measurement between two points. The GEOALIGN® Marking System may also facilitate geographic alignment of an adjunctive therapy that includes the same GEOALIGN® Marking System.

**NOTE:** The GEOALIGN® Marking System provides an approximation that may not be an exact representation of the actual distance traveled intravascularly and should be confirmed under fluoroscopy.


Figure 2

**NOTE:** The GEOALIGN® Marking System includes non-radiopaque white markings and are designed to be utilized outside the sheath.

#### Indications for Use

The ULTRAScore™ Focused Force PTA Balloon is intended to dilate stenoses in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also recommended for post dilatation of balloon expandable stents, self-expanding stents, and stent grafts in the peripheral vasculature.

#### Contraindications

The ULTRAScore™ Focused Force PTA Balloon is contraindicated:

- Where there is the inability to cross the target lesion with a guidewire
- For use in the coronary or neuro vasculature

#### Warnings

1. Contents supplied STERILE using ethylene oxide (EO). Non-Pyrogenic. Do not use if sterile barrier is opened or damaged. Do not reuse, reprocess or re-sterilize. Use the catheter prior to the "Use By" date specified on the package label.
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminate period of time. The residue of biological material can promote the contamination of the device with

pyrogens or microorganisms which may lead to infectious complications.

3. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminate degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. To reduce the potential for vessel damage or difficulty in deflating, the inflated diameter and length of the balloon should approximate the diameter and length of the vessel just proximal and distal to the stenosis.
5. When the catheter is exposed to the vascular system, the location of the balloon should be confirmed while under high quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip or catheter breakage, catheter kink, or balloon separation.
6. Do not exceed the RBP recommended for this device. Balloon rupture or difficulty in deflation may occur if the RBP rating is exceeded. To prevent over pressurization, use of a pressure monitoring device is recommended.
7. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

#### Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. The ULTRAScore™ Focused Force PTA Balloon should only be used by physicians experienced in the performance of percutaneous transluminal angioplasty.
3. It is recommended to consider the use of anti-coagulants, anti-platelet agents, and/or vasodilators in conformance with the accepted standard of practice or institutional guidelines surrounding peripheral endovascular procedures.
4. For ULTRAScore™ .014" guidewire sizes only, in order to activate the hydrophilic coating, wet the ULTRAScore™ balloon and catheter with sterile saline or wipe the balloon catheter with sterile saline saturated gauze immediately prior to its insertion in the body. Do not wipe the balloon catheter with dry gauze.
5. The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the PTA catheter through a smaller size introducer sheath than indicated on the label.
6. Use the recommended balloon inflation medium (25% contrast medium/75% sterile saline solution). Never use air or other gaseous medium to inflate the balloon.
7. The ULTRAScore™ Focused Force PTA Balloon should be used with caution for procedures involving calcified lesions, stents or synthetic vascular grafts due to the abrasive nature of these lesions.
8. Fully evacuate the balloon prior to withdrawing the system. Larger sizes of ULTRAScore™ Focused Force PTA Balloons may exhibit slower deflation times.
9. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath, determine if contrast medium is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the sheath and then completely evacuate the contrast before proceeding to withdraw the balloon.
10. If resistance is still felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and guidewire/introducer sheath as a single unit.
11. Do not continue to use the balloon catheter if the shaft has been bent or kinked.
12. For ULTRAScore™ .014" guidewire sizes only, prior to re-insertion through the introducer sheath, re-activate the hydrophilic coating, and clean the balloon catheter by wiping the balloon catheter with sterile saline saturated gauze and rinsing with sterile saline. Do not wipe the balloon catheter with dry gauze.
13. GEOALIGN® Marking System is designed to be used as an additional reference tool to accompany the interventionalist standard operation procedure. The use of fluoroscopic imaging is recommended following positioning of the catheter to the target lesion and prior to balloon deployment.

#### Potential Adverse Reactions

The complications that may result from a peripheral balloon dilatation procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stroke

- Thrombosis
- Vessel dissection, perforation, rupture, or spasm

## Directions for Use

### Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultraviolet light sources. Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date. Do not use if packaging is damaged or opened.

### Equipment Required

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (10 ml or larger)
- Appropriate introducer sheath and dilator set
- .014" or .035" guidewire
- Optional: Catheter stabilization device

### Dilatation Catheter Preparation

1. Remove the balloon catheter from package. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
2. Remove the stylet and protective sheath by grasping the balloon catheter just proximal to the balloon and with the other hand, gently grasp the mandrel and sheath protector and slide distally off of the balloon catheter.

**Precaution:** While removing the protective sheath, care should be taken to avoid damaging the balloon catheter.

3. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium (25% contrast medium/75% sterile saline solution). Do not use air or any gaseous medium to inflate the balloon.
4. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
5. Connect the syringe to the stopcock.
6. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
7. Repeat step #6 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
8. In order to activate the coating, wet the balloon catheter with sterile saline or wipe the balloon catheter with saline saturated gauze immediately prior to its insertion in the introducer sheath.
9. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.

**Precaution:** Do not wipe the balloon catheter with dry gauze.

### Use of the ULTRAScore™ Focused Force PTA Balloon

1. Backload the distal tip of the ULTRAScore™ Focused Force PTA Balloon over the pre-positioned guidewire and advance the tip to the introduction site.
2. Advance the catheter through the introducer sheath/guide catheter and over the wire to the site of inflation. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.

**Note:** If using a hydrophilic guidewire, ensure that it is kept hydrated with sterile normal saline at all times.

3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place and slowly inflate the balloon with an inflation device.

**Note:** Balloon lengths 100mm and longer have two marker bands located at the distal portion of the working length of the balloon and one marker band at the proximal portion. The working length of the balloon is measured at the proximal edge of the proximal marker band to the distal edge of the most distal marker band.

**Precaution:** Do not inflate the balloon past its labeled rated burst pressure.

4. Apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy and that no contrast is left in the balloon.
5. While maintaining negative pressure and the position of the guidewire, grasp the balloon catheter just outside the introducer sheath/guide catheter and withdraw the deflated dilatation catheter over the wire through the introducer sheath/guide catheter.

**Optional:** To ensure that the sheath does not move during the procedure, use a catheter stabilization device to stabilize the sheath at the access site.

**Optional:** When used as a location reference tool, utilize the GeoALIGN® Marking System to assist with measuring how far the catheter has advanced. The GeoALIGN® Marking System provides an approximation that may not be an exact representation of the actual distance traveled intravascularly and should be confirmed under fluoroscopy (Reference Figure 3).


Figure 3

**Optional:** When used for intravascular length measurement, advance or retract the catheter between Points A and B utilizing the GeoALIGN® Marking System to measure the difference. The GeoALIGN® Marking System provides an approximation that may not be an exact representation of the actual distance traveled intravascularly and should be confirmed under fluoroscopy (Reference Figure 4).


Figure 4

**Optional:** To facilitate repeat catheter alignment with an additional device with the GeoALIGN® Marking System, ensure the GeoALIGN® markings are at the same location outside the sheath as the initial device.

**Note:** Ensure the sheath does not move during the withdrawal of the initial device or the advancement of subsequent devices.

### Balloon Reinsertion

**Precaution:** Do not continue to use the balloon catheter if the shaft has been bent or kinked.

**Precaution:** For ULTRAScore™ .014" guidewire sizes only, prior to re-insertion through the introducer sheath, re-activate the hydrophilic coating, and clean the balloon catheter by wiping the balloon catheter with sterile saline saturated gauze and rinsing with sterile saline. Do not wipe the balloon catheter with dry gauze.

1. Load the balloon catheter onto a guidewire.
2. Advance the balloon catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon catheter.
3. Continue the procedure according to the "Use of the ULTRAScore™ Focused Force PTA Balloon" section herein.

**Warning:** After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

**Warning:** After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

### Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, at Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product is not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**Copyright © 2018 C. R. Bard, Inc. All rights reserved.**
Copyright © 2018 C. R. Bard, Inc. Tous droits r serv s.
Copyright © 2018 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2018 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2018 C. R. Bard, Inc. Todos los derechos reservados.
Copyright © 2018 C. R. Bard, Inc. Alle rechten voorbehouden.