# Modified ONFLEX™ Mesh

Self-Expanding Lightweight Mesh with Strap  
for Open Preperitoneal Inguinal Hernia Repair

## Instructions for Use

Mode d'emploi  
Gebrauchsanweisung  
Istruzioni per l'uso  
Instrucciones de uso  
Gebruiksaanwijzing  
Instruções de utilização  
Οδηγίες Χρήσης  
Brugervejledning  
Bruksanvisning  
Käyttöohjeet  
Bruksanvisning  
Instrukcja użycia  
Használati útmutató  
Návod k použití  
Kullanım Talimatları  
Инструкции по применению  
使用説明  
사용 설명서

 Single Use

Partially Absorbable

R<sub>x</sub> only

STERILE EO

Manufacturer:

Davol Inc.  
Subsidiary of C. R. Bard, Inc.  
100 Crossings Boulevard  
Warwick, RI 02886 USA  
1-401-825-8300 • 1-800-556-6275

EC REP

  
0086

Bard Limited  
Crawley, UK  
RH11 9BP

Medical Services & Support  
Clinical Information Line  
1-800-562-0027

DAVOL INC.

PK3798862 153R

# Modified ONFLEX™ Mesh

Self-Expanding Lightweight Mesh with Strap  
for Open Preperitoneal Inguinal Hernia Repair

## DESCRIPTION:

The Modified ONFLEX™ Mesh is a self-expanding, non-absorbable, sterile prosthesis made from monofilament polypropylene mesh and has a lightweight large pore design. This construction allows a prompt fibroblastic response through the interstices of the mesh as observed in a preclinical model, which may not correlate to performance in humans. The Modified ONFLEX™ Mesh has an anatomical shape designed to cover potential defect areas.

The Modified ONFLEX™ Mesh also contains two pockets to facilitate insertion and positioning of the device. The positioning pockets are located on the larger medial apex of the mesh and the lateral apex of the mesh. In addition to the pocket, the mesh also contains straps to facilitate positioning and fixation of the device. The Modified ONFLEX™ Mesh comes packaged with an onlay which is available in one size and is optional based on surgeon preference.

The device contains SORBAFLEX™ Memory Technology which provides memory and stability to the device, facilitating ease of initial insertion and proper placement of the device. The SORBAFLEX™ Memory Technology is comprised of an extruded polydioxanone (PDO) absorbable monofilament that is contained within a knitted polypropylene mesh tube. The PDO monofilament fully degrades *in vivo* by means of hydrolysis and has been found to elicit an inflammatory response during absorption. Absorption is essentially complete in 6–8 months. The PDO monofilament is dyed violet by adding D & C Violet No.2.

## INDICATIONS:

The Modified ONFLEX™ Mesh is indicated for use in the reinforcement of soft tissue where weakness exists, such as in the repair of inguinal hernias.

## CONTRAINDICATIONS:

- Use of this device is contraindicated for infants, children, or pregnant women, whereby future growth will be compromised by use of such mesh material.
- Literature reports that there is a possibility for adhesion formation when polypropylene is placed in direct contact with the bowel or viscera.

## WARNINGS:

- The use of any synthetic mesh or patch in a contaminated or infected wound can lead to fistula formation and/or extrusion of the mesh and is not recommended.
- If an infection develops, treat the infection aggressively. Consideration should be given regarding the need to remove the mesh. Unresolved infection may require removal of the mesh.
- If unused mesh has been in contact with instruments or supplies used on a patient or contaminated with body fluids, discard mesh with care to prevent risk of transmission of viral infections.
- To prevent recurrences when repairing hernias, the mesh should be large enough to provide sufficient overlap beyond the margins of the defect. If fixation is used, careful attention to mesh fixation placement and spacing will help prevent excessive tension or gap formation between the mesh and fascial tissue.
- To avoid injury, careful attention is required if fixating the mesh in the presence of nerves, vessels, or spermatic cord structures.
- The mesh is supplied sterile. Inspect the packaging to be sure it is intact and undamaged prior to use.
- This device is for single use only. Do not resterilize or reuse any portion of the Modified ONFLEX™ Mesh.
- Do not cut or reshape the Modified ONFLEX™ Mesh, except for the positioning strap, as this could affect its effectiveness. Care should be taken not to cut or nick the SORBAFLEX™ PDO monofilament.

- Excess positioning strap material above the fixation point must be cut off at the level of the fascia and discarded to eliminate excess material from remaining in the body.

## PRECAUTIONS:

- Please read all instructions prior to use.
- Only physicians qualified in the appropriate surgical techniques should use this prosthesis.
-  Care should be taken not to cut or nick the SORBAFLEX™ PDO monofilament.
- If fixation is used, care should be taken to ensure that the mesh is adequately fixated. If necessary, additional fasteners and/or sutures should be used.

## ADVERSE REACTIONS:

Possible complications may include, but are not limited to, seroma, adhesion, hematoma, pain, infection, inflammation, extrusion, erosion, migration, fistula formation and recurrence of the hernia or soft tissue defect. If the SORBAFLEX™ PDO monofilament is cut or damaged, additional complications may include, but are not limited to, bowel or skin perforation and infection.

## INSTRUCTIONS FOR USE:

Select the appropriate size of the Modified ONFLEX™ Mesh. Once the hernia has been reduced, dissect the appropriate anatomical space to accommodate the selected size of the Modified ONFLEX™ Mesh. Using the pockets, roll the mesh without folding the PDO monofilament along the long axis and guide the mesh into the previously created space while maintaining the free end of the positioning strap outside of the body. The positioning strap is manipulated to facilitate proper placement and positioning of the mesh. Gently pulling up on the positioning strap will allow the mesh to flatten against the inguinal floor. To secure the mesh, the mesh positioning strap is pulled up and the straps are secured to the margins of the defect. Excess positioning strap material above the fixation points must be cut off at the level of the fascia and discarded. When more mesh coverage is desired, the optional onlay may be used for additional reinforcement. Close wound.

## TRACEABILITY:

A traceability label which identifies the type, size, expiration date, and lot number of the prosthesis is attached to every package. This label should be affixed to the patient's permanent medical record to clearly identify the device which was implanted.

If you experience a product failure, please contact Davol, Inc. at 1-800-556-6275 for instructions on returning the product.

Bard, Davol, Onflex and SorbaFlex are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 2014 C. R. Bard, Inc. All Rights Reserved.

	Medium
	Large
	Contents
	Do not use if package is damaged.
	Actual Size
	Partially Absorbable
	U. S. Federal law restricts this device to sale by or on the order of a physician
	Peel Here

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Copyright © 2014 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2014 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2014 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2014 C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2014 C. R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2014 C. R. Bard, Inc. Todos os direitos reservados.