# ONFLEX™ Mesh

Self-Expanding Lightweight Mesh

for Open Preperitoneal Inguinal Hernia Repair

---

## Instructions for Use

---

Mode d'emploi  
Gebrauchsanweisung  
Istruzioni per l'uso  
Instrucciones de uso  
Gebruiksaanwijzing  
Instruções de utilização  
Οδηγίες Χρήσης  
Brugervejledning  
Bruksanvisning  
Käyttöohjeet  
Bruksanvisning  
Instrukcja użycia  
Használati útmutató  
Návod k použití  
Kullanım Talimatları  
Инструкции по применению  
使用説明  
사용 설명서

---

 Single Use

Partially Absorbable

R<sub>x</sub> only

**STERILE EO**

---

Manufacturer:

Davol Inc.  
Subsidiary of C. R. Bard, Inc.  
100 Crossings Boulevard  
Warwick, RI 02886 USA  
1-401-825-8300 • 1-800-556-6275

**EC REP**

Bard Limited  
Crawley, UK  
RH11 9BP

**CE**  
0086

Medical Services & Support  
Clinical Information Line  
1-800-562-0027

**BAIRD**

DAVOL INC.

PK3798861 153R

# ONFLEX™ Mesh

Self-Expanding Lightweight Mesh  
for Open Preperitoneal Inguinal Hernia Repair

## DESCRIPTION:

The ONFLEX™ Mesh is a self-expanding, non-absorbable, sterile prosthesis made from monofilament polypropylene mesh and has a lightweight large pore design. This construction allows a prompt fibroblastic response through the interstices of the mesh as observed in a preclinical model, which may not correlate to performance in humans. The ONFLEX™ Mesh has an anatomical shape designed to cover potential defect areas. The ONFLEX™ Mesh also contains a pocket on the larger medial apex of the mesh to facilitate insertion and positioning of the device.

The device contains SORBAFLEX™ Memory Technology which provides memory and stability to the device, facilitating ease of initial insertion and proper placement of the device. The SORBAFLEX™ Memory Technology is comprised of an extruded polydioxanone (PDO) absorbable monofilament that is contained within a knitted polypropylene mesh tube. The PDO monofilament fully degrades *in vivo* by means of hydrolysis and has been found to elicit an inflammatory response during absorption. Absorption is essentially complete in 6-8 months. The PDO monofilament is dyed violet by adding D & C Violet No.2.

## INDICATIONS:


The ONFLEX™ Mesh is indicated for use in the reinforcement of soft tissue where weakness exists, such as in the repair of inguinal hernias.

## CONTRAINDICATIONS:

1. Use of this device is contraindicated for infants, children, or pregnant women, whereby future growth will be compromised by use of such mesh material.
2. Literature reports that there is a possibility for adhesion formation when polypropylene is placed in direct contact with the bowel or viscera.

## WARNINGS:

1. The use of any synthetic mesh or patch in a contaminated or infected wound can lead to fistula formation and/or extrusion of the mesh and is not recommended.
2. If an infection develops, treat the infection aggressively. Consideration should be given regarding the need to remove the mesh. Unresolved infection may require removal of the mesh.
3. If unused mesh has been in contact with instruments or supplies used on a patient or contaminated with body fluids, discard mesh with care to prevent risk of transmission of viral infections.
4. To prevent recurrences when repairing hernias, the mesh should be large enough to provide sufficient overlap beyond the margins of the defect. If fixation is used, careful attention to mesh fixation placement and spacing will help prevent excessive tension or gap formation between the mesh and fascial tissue.
5. To avoid injury, careful attention is required if fixing the mesh in the presence of nerves, vessels, or spermatic cord structures.
6. The mesh is supplied sterile. Inspect the packaging to be sure it is intact and undamaged prior to use.
7. This device is for single use only. Do not resterilize or reuse any portion of the ONFLEX™ Mesh.
8.  Do not cut or reshape the ONFLEX™ Mesh, except at the opening in the interrupted SORBAFLEX™ PDO monofilament, to accommodate the spermatic cord and outside the blue limit line in the lateral portion of the mesh, as this could affect its effectiveness. Care should be taken not to cut or nick the SORBAFLEX™ PDO monofilament.

## PRECAUTIONS:

1. Please read all instructions prior to use.
2. Only physicians qualified in the appropriate surgical techniques should use this prosthesis.
3.  Care should be taken not to cut or nick the SORBAFLEX™ PDO monofilament.
4. If fixation is used, care should be taken to ensure that the mesh is adequately fixated. If necessary, additional fasteners and/or sutures should be used.

## ADVERSE REACTIONS:

Possible complications may include, but are not limited to, seroma, adhesion, hematoma, pain, infection, inflammation, extrusion, erosion, migration, fistula formation and recurrence of the hernia or soft tissue defect. If the SORBAFLEX™ PDO monofilament is cut or damaged, additional complications may include, but are not limited to, bowel or skin perforation and infection.

## INSTRUCTIONS FOR USE:

Select the appropriate size of the ONFLEX™ Mesh. Once the hernia has been reduced, dissect the preperitoneal space to accommodate the selected size of the mesh. Using the pocket, roll the mesh along the long axis without folding the PDO monofilament and guide the medial apex of the mesh into the preperitoneal space to adequately cover the defect. The lateral apex of the mesh can be placed in either the preperitoneal space or onlay position above the muscle layer to adequately cover the defect. Based on surgeon preference and surgical technique, a keyhole or slit may be cut into the mesh at the area of the interrupted PDO monofilament and the tails of the mesh may be placed around the spermatic cord. The tails can be brought together with sutures, while ensuring adequate space to accommodate the spermatic cord. Secure the mesh and close wound.


## TRACEABILITY:

A traceability label which identifies the type, size, expiration date, and lot number of the prosthesis is attached to every package. This label should be affixed to the patient's permanent medical record to clearly identify the device which was implanted.

If you experience a product failure, please contact Davol, Inc. at 1-800-556-6275 for instructions on returning the product.

Bard, Davol, Onflex and SorbaFlex are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 2014 C. R. Bard, Inc. All Rights Reserved.

	Medium
	Large
	Contents
	Do not use if package is damaged.
	Actual Size
	Partially Absorbable
	U. S. Federal law restricts this device to sale by or on the order of a physician
	Peel Here

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Copyright © 2014 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2014 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2014 C. R. Bard, Inc. Tutti i diritti riservati.
Bard, Davol, Onflex y SorbaFlex son marcas comerciales y/o marcas comerciales registradas de C. R. Bard, Inc. o de una filial. Copyright © 2014 C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2014 C. R. Bard, Inc. Alle rechten voorbehouden.
Bard, Davol, Onflex e SorbaFlex são marcas comerciais e/ou marcas registadas da C. R. Bard, Inc. ou de uma empresa afiliada. Copyright © 2014 C. R. Bard, Inc. Todos os direitos reservados.