# XENMATRIX™ Surgical Graft

Porcine Dermal Matrix

---

## INSTRUCTIONS FOR USE

---

Mode d'emploi  
Gebrauchsanweisung  
Istruzioni per l'Uso  
Instrucciones de uso  
Gebruiksaanwijzing  
Instruções de Utilização  
Οδηγίες χρήσης  
Brugervejledning  
Bruksanvisning  
Käyttöohjeet  
Bruksanvisning  
Instrukcja użycia  
Használati útmutató  
Návod k použití  
Kullanım Talimatları  
使用説明  
사용 설명서  
Инструкции по применению

Do Not Resterilize

Single Use

R<sub>x</sub> Only

STERILE R

**BAIRD**

DAVOL INC.

# XENMATRIX™ Surgical Graft

Porcine Dermal Matrix

Please read all instructions prior to use.

## PRODUCT DESCRIPTION

XENMATRIX™ Surgical Graft is an acellular, sterile, non-pyrogenic porcine dermal matrix packed hydrated in sterile saline for use in the reconstruction of soft tissue deficiencies. The XENMATRIX™ Surgical Graft is an AQUAPURE™ processed material.

## STORAGE

Store the XENMATRIX™ Surgical Graft in a dry environment and protected from direct sunlight. Store refrigerated or at a temperature not to exceed 25°C (77°F). It is for single use only. Do not use if the package is damaged or open, or if the center of the temperature indicator on the foil pouch is black.

## INDICATIONS

Intended for implantation to reinforce soft tissue where weakness exists and for surgical repair of damaged or ruptured soft tissue, including: abdominal plastic and reconstructive surgery; muscle flap reinforcement; hernia repair including abdominal, inguinal, femoral, diaphragmatic, scrotal, umbilical, and incisional hernias.

## CONTRAINDICATIONS

1. XENMATRIX™ Surgical Graft should not be used on patients with known sensitivity to porcine products.
2. Not for reconstruction of cardiovascular defects.
3. Not for reconstruction of central nervous system or peripheral nervous system defects.
4. Use of this product in applications other than those indicated has the potential for serious complications.

## WARNINGS

1. **This device has been designed for single use only. Reuse, reprocessing, resterilization or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device and may lead to device failure which may result in injury to the patient.**

Reuse, reprocessing, resterilization or repackaging may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient or end user.

2. **DO NOT RESTERILIZE.**
3. **Prior to use, carefully examine package and product to verify neither is damaged and that all seals are intact. Do not use if the package is damaged or open, or if the center of the temperature indicator on the foil pouch is black.**
4. **If an infection develops, it should be treated aggressively.**
5. **An allergic reaction, which is unrelated to other therapy, is an indication to consider removal of XENMATRIX™ Surgical Graft.**
6.  **After use, any unused product and packaging should be treated as a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.**

## PRECAUTIONS

1. Please read all instructions prior to use.
2. Only physicians qualified in the appropriate surgical techniques should use this surgical graft.
3. Strict aseptic technique should be followed.
4. U. S. federal law restricts this device to sale by or on the order of a physician.
5. The surgeon should thoroughly understand the surgical procedure and the performance characteristics of the surgical graft.
6. Place device in maximum possible contact with healthy, well-vascularized tissue to promote cell ingrowth and tissue remodeling.
7. When unable to close skin over the XENMATRIX™ Surgical Graft, ensure that the implant remains moist. Avoid drying of the implant through "continued suction devices" as this may negatively impact the performance of the implant.

## ADVERSE REACTIONS

Potential complications with the use of any prosthesis may include, but are not limited to, allergy, seroma, infection, inflammation, adhesion, fistula formation, hematoma and recurrence of tissue defect.

## INSTRUCTIONS FOR USE

Use sterile gloves and/or atraumatic instruments to open package.

### Preparation

XENMATRIX™ Surgical Graft is ready for use out of the package. The graft may be placed in sterile saline if desired prior to implantation.

### Sizing

Cutting XENMATRIX™ Surgical Graft to the proper size is essential. The surgical graft should be large enough to extend beyond the margins of the defect. A minimum overlap of 3–5 cm is recommended. Cutting XENMATRIX™ Surgical Graft too small for the defect can cause excessive tension on the suture line.

### Deployment

XENMATRIX™ Surgical Graft may inhibit deployment through trocars during laparoscopic repairs. If the XENMATRIX™ Surgical Graft cannot be easily deployed down the trocar, remove the trocar and insert the implant through the incision. Reinsert trocar.

### Fixation

XENMATRIX™ Surgical Graft should be fixated with suture under minimal or no tension. Placement of sutures 2–3 cm from the edge of the graft is recommended.

## TRACEABILITY

A traceability label that identifies the type, size, and lot number of the implant is attached to every package. The label should be incorporated into the patient's permanent medical record to clearly identify the device that was implanted.

AquaPure, Bard, Davol and XenMatrix are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 2009-2014 C. R. Bard, Inc. All Rights Reserved

	Rectangle		Contents
	Square		Do not use if package is damaged.
	U. S. federal law restricts this device to sale by or on the order of a physician.		

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Copyright © 2009-2014 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2009-2014 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2009-2014 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2009-2014 C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2009-2014 C. R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2009-2014 C. R. Bard, Inc. Todos os direitos reservados.