

PRESCRIPTIVE INFORMATION

INDICATIONS FOR USE

The SUREPOINT® (PT STEPPER™) equipment is a precision stepping & stabilization system designed to securely hold the designated ultrasound transducer. It allows precise manual rotation and stepwise movement of the transducer along the longitudinal axis. It is available in a portable, table and floor configuration.

CAUTIONS

- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- DO NOT AUTOCLAVE. If placed in an autoclave, the SUREPOINT® (PT STEPPER™) equipment will be severely damaged and the warranty will be null and void.