

PRESCRIPTIVE INFORMATION

INDICATIONS

- For placement of radioactive seeds in transperineal implant procedures.

CAUTIONS

- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- BARD® BRACHYSTAR® Applicator Seed and FASTFILL® Implant Needle - Using the stylet with excessive force to manipulate lodged seeds may damage the seed.
- BARD® BRACHYSTAR® FASTFILL® Implant Needle - Needles are not intended to penetrate bone. If resistance is encountered, verify needle position. Do not push into bone; this may cause needle to bend or break. Replace needle if cannula, or point, is damaged.

WARNINGS

- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and with applicable laws and regulations.

NOTES

- Sterile unless package has been opened or damaged.
- Do not use if package is opened or damaged.