

# PRESCRIPTIVE INFORMATION

**INDICATIONS FOR USE:** For catheter irrigation.

## **CAUTIONS**

- Secure tip tightly prior to use.
- Avoid excessive syringe pressure when using luer adapter, as adapter may dislodge.

## **INFORMATION**

- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.
- Sterile unless package is opened or damaged.
- Do not use if package is opened or damaged.
- Single use.