



## PRESCRIPTIVE INFORMATION

Use: For urological use only.

### Caution

- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- This product contains natural rubber latex which may cause allergic reactions.

### Warnings

- Do not use ointments or lubricants having a petroleum base. This will damage latex. Visually inspect the product for any imperfections or surface deterioration prior to use.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

### Notes

- Single patient use, Do not reuse.
- Sterile unless package is opened or damaged.

For the latest information, always check the "Instructions for Use" that comes packaged with the product. This site is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider.

### TRAINING AND EDUCATION

For additional training and education on a number of topics, visit our Training Center

[read more ▶](#)

### CONTACT US

[read more ▶](#)

[give site feedback](#) 

[Home](#) [Terms of Use](#) [Privacy Policy](#) [Contact Us](#) [C.R. Bard](#)

Copyright © 2013 C. R. Bard, Inc. All rights reserved. All products on this website may not be available for sale in all geographies. Please contact your local Bard Medical representative for information about products available in your area.

This site is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider.