

HemoSplit*
LONG-TERM HEMODIALYSIS CATHETER

HemoSplit*^{XX}
LONG-TERM HEMODIALYSIS CATHETER

with



BIOBLOC*
SILVER SULFADIAZINE COATING

If you could reduce bacterial adhesion to hemodialysis catheters, why wouldn't you?



BAIRD
ACCESS SYSTEMS

*Advancing the Delivery of Health Care.**

The weapon against bacteria so you have one less worry.



It goes without saying— Bacterial adhesion is the last thing you want to battle when delivering patient care. Catheter failure due to bacterial adhesion may disrupt your workflow, create high costs for your facility, and may create more pain and discomfort to your patient.


Now you can be confident in your ability to reduce bacterial adhesion to the catheter, as supported by *in-vitro* testing, when you use the **HemoSplit*** and **HemoSplit* XK** Catheters with **BioBloc** coating from Bard Access Systems. **BioBloc** coating is innovative technology that uses silver sulfadiazine that has a proven history in the medical device industry, most often recognized in burn creams.

As well, **BioBloc** coating is—

- **Simple to use**
- **Safe for the patient**
- **Affordable enough to use on virtually all patients**

21 Days of Protection —

Bard Access Systems has developed **BioBloc** coating that reduces bacterial adhesion to the catheter by 99.9% in the catheter tunnel for a period of 21 days as tested in an *in-vitro* model against *Staphylococcus epidermidis*, *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Candida albicans*, *Enterococcus faecalis* and *Escherichia coli*. And because silver sulfadiazine is not an antibiotic, it is rare for bacteria to develop resistance.



- Ordering Information -

16 Fr. HemoSplit XK Catheter with BioBloc Coating	14.5 Fr. HemoSplit Catheter with BioBloc Coating				
	[568] - Standard Kit	[573] - Standard Kit	[574] - Microintroducer Kit	Tip to Cuff	Tip to Bifurcation
Straight	Straight				
	5733153			15 cm	20 cm
5683693	5733693		5743693	19 cm	24 cm
5683733	5733733		5743733	23 cm	28 cm
5683273	5733273		5743273	27 cm	32 cm
5683313	5733313		5743313	31 cm	36 cm
5683353	5733353			35 cm	40 cm
5684423	5734423			42 cm	47 cm

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

Learn More -

To learn more about **HemoSplit** and **HemoSplit XK** Catheters with **BioBloc** coating, contact your Bard Access Systems representative, call us toll-free at **1-800-545-0890**, or visit us on the web at **www.bardaccess.com**.



ACCESS SYSTEMS

Bard Access Systems, Inc.

Salt Lake City, Utah 84116 USA
801-595-0700
www.bardaccess.com

Customer Service: 800-545-0890
Clinical Information: 800-443-3385